




MEDICAL SECTOR

INFATN BILIRUBIN PHOTOTHERAPY



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Infant Bilirubin Phototherapy
Infant Bilirubin Phototherapy Unit
Model YR02194



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Infant Bilirubin Phototherapy Infant Bilirubin Phototherapy Unit Model YR02194



Product Description:

The newborn patients with hyperbilirubinemia all need best care for threatening jaundice. Our line of phototherapy products supply you the reliability and stability intensive phototherapy by LED light for your baby patient. The unit has a long-life that provides 20,000 to 0,000 hours of light therapy.

Features:

- ✓ LED fluorescent lamp as the phototherapy radiating elements
- ✓ Angle of lamp head can be adjusted horizontally.
- ✓ Height of lamp head can be adjusted.
- ✓ Stainless steel support stand can be adjusted.
- ✓ Wheels can be locked
- ✓ Timer to record the phototherapy operation.
- ✓ Use in incubator, irradiant heat bassinets, and any kind of bassinet.

Application:

The YR02194 unit like YR02194 unit with blue florescent tube has independent framework can use with incubator, irradiant heat bassinets, and any kind of bassinet in family care. Widely applicable in all hospitals and clinics and optimize for jaundiced baby..



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
 KALSTEIN FRANCE, S. A. S



KALSTEIN® Infant Bilirubin Phototherapy
Infant Bilirubin Phototherapy Unit
Model YR02194



Technical Description:

Model	YR02194
Irradiation source	LED Light
Irradiation wavelength range (300 ~ 550) nm	•
Irradiance adjustable from 1 to 99 levels	3 000 uw / cm2
Time calculator	•
Countdown	•
Total time	•
Display range 0 – 9999 h	•
Adjustable working height	•
Effective shelf life	5 000h
Nominal Lifespan of Blue Light	20 000h
Energy consumption	≤ 50W
AC motor	•
Function	
Irradiation wavelength range	400 ~ 550nm
Adjustable irradiation intensity (max.)	3 000 uw / cm ²
Wavelength of wave crest	460nm
Irradiance adjustable from 1 to 99 soft	•
Hourly	•
A fixed time	•
Time totalizer	•
Display range (hour)	0 – 9 999h
Structural parts	
Adjustable in vertical height	1 100 – 1 589mm
Useful life	5 000h
Nominal Lifespan of Blue Light	20 000h
Dimension	500 x 600 x (1 355 ~ 1 705) mm
Electric parts	
Power supply	AC (220 ± 22) V (50 ± 1) Hz
Energy consumption	≤50W
Room temperature	10 – 40 ° C
Relative humidity	≤80%
Motor	AC



#Letsgivemore ♥
 Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

