


# Kalstein A different accompaniment, at your service

# **LABORATORY LINES**

SAFETY STORE CABINETS





All rights reserved ® KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,

• 34000 Montpellier, France.







OUR SERVICES

### **Benefits and Support**

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



# Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



#### **Quick Response**

Our work team is always available to response all your consults or questions, in order to support you in any situation.





#### #Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



#### **Technical Support**

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.





#### **Delivery Logistics**

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or.air.



#### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



#### **90** Gallon Explosion-proof Cabinet

#### Model YR05472//YR05472-2





#### **Product Description**

- ✓ Environmental protection spray, strong insulation
- ✓ Three-point linkage, double lock
- ✓ Strong load-bearing, adjustable laminate
- ✓ Electrostatic grounding, safety upgrade
- ✓ Three colors, each with its own strengths

#### **Technical Description**

Model	YR05472	YR05472-1	YR05472-2
Dimensions (H x W x D)	1 650 x 1 090 x 865		
Weight (Kg)	165		
Volume	90 GAL / 340 GAL		
Door Type	Double door / Manual		
Description	Yellow for Flammable Liquids	Red for Corrosive liquids	Blue for Combustible liquids
Adjustable shelf	2 Layer		



#### #Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.





34000 Montpellier, France.
 Tit: +33 467158849 / +33 680760710/ +33 663810023
 https://kaidein.eu

https://kaistein.eu Kalstein France, S.A. S



#### **90** Gallon Explosion-proof Cabinet

#### Model YR05472//YR05472-2




#### **Features**

•Double lock safe and reliable: Two locks, two managers. Double security, pad lock and cabinet lock, stardarized management, improve safety factor



 Three-point linkage solid lock post: Three-point linkage cabinet door lock tongue locking structure, upper, middle, and lower three-point linkage solid core lock columns to increase the security
 of
 the
 cabinet



 Adjustable fire vent: The fire vents on both sides of the cabinet body are equipped with adjustable valves, which can rotate to control the air volume and reduce the risk of harmful gas accumulation in the cabinet







#### #Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved @ KALSTEIN France S. A. S., Optimum Business Center 450 Rue Baden Powell,

• 34000 Montpellier, France.

Tif: +33 467158849 / +33 680760710/ +33 663810023

https://kaistein.eu Kalstein France, S.A.S



#### **90** Gallon Explosion-proof Cabinet

#### Model YR05472//YR05472-2



 Liquid leakproof galvanized shelves: 4-110 gallons, regardless of size, the shelves are designed with a leaked liquid guide groove, which is not afraid of leakage. The laminate can be adjusted up and down according to the demand





 Anti-static: each model is equipped with an OSHA standard electrostatic prevention wire clip. Ensure safe use and management






#### #Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.





• 34000 Montpellier, France. Tif: +33 467158849 / +33 680760710/ +33 663810023

https://kaistein.eu Kalstein France, S.A. S







#LETSGIVEMORE >

# WITH THE ADCQUISITION OF A KALSTEIN EQUIPMENT

YOU MAKE A CONTRIBUTION TO:





https://onetreeplanted.org



Solidarité et coopération biomédicale

https://www.humatem.org/



https://www.maniapurefoundation.org/



# ANY DOUBT? Contact us!

#### **Articles Kalstein**

(Social Networks and Official Websites)







Article of the week: How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites





**Article of the week:** What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites https://kalstein.eu/articles

#### **Kalstein Training Videos**

(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



### MONTPELLIER - FRANCE HEADQUARTERS

Optimum Business Center, 450 Rue Baden Powell, 34000 Montpellier - France.

Fax: +33 (0) 4 67 15 88 49 Tlf: +33 (0) 6 80 76 07 10 sales@kalstein.eu

https://kalstein.net/ https://www.kalstein.fr/

#### Official Websites for Latin America.

