

LABORATORY SECTOR

REFRACTOMETER







All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S





OUR SERVICES Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings In any part of the world, receive your

induction or training from our specialized team of engineers



Quick Response Our work team is always available to response all your consults or questions, in order to support you in any situation.





#Letsgivemore 💙

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or.air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.





KALSTEIN® Refractometer ATC Wine Refractometers Model YR05931 // YR05936





#Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.





All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S



www.kalstein.eu

KALSTEIN® Refractometer ATC Wine Refractometers Model YR05931 // YR05936



Product Description:

The Wine Refractometer is a metal device that provides you with fast & accurate results. It's dual-scale provides you the direct reading of the alcohol as well as brix measurements. FULLY EQUIPPED with ATC that automatically adjust itself to correct temperature discrepancies during use. Refractive index is very temperature dependent; it is important to use a Refractometer with automatic temperature compensation to give you the most accurate measurement in any type of environment temperature. Capable of measuring with high accuracy and provides precise results and is equipped with vivid and sharp reticle chart for easy and comfortable reading.

Features:

- ✓ High quality body construction
- ✓ Uses ambient light only which means battery or power source is not required
- ✓ Easy to focus and calibrate
- ✓ Cushioned with soft & comfortable non-slip rubber
- ✓ With ATC function(Auto Temperature Compensation)
- ✓ Durable and built to last long

Technical Specifications:



Model	Description	Measurement Range	Min. Div.	Accuracy	ATC	ATC Temp.	Material
YR05931	ATC alcohol 0~80% v/v optical refractometer	0-80%v/v	1%v/v	±1%v/v	NO	_	Pure aluminium
YR05932	ATC Wine 0-25% Vol 0-40%Brix optical refractometer	0-25%Vol 0- 40%Brix	0.2%Vol 0.2%Brix	±0.2%Vol ±0.2%Brix	Yes	10-30° C	Pure aluminium
YR05933	0-25%Vol 0-170 °Oe	0-25%Vol 0-170 °Oe	0.2%Vol 1°Oe	±0.2%Vol ±1°0e	Yes	10-30° C	Pure aluminium
YR05934	ATC Wine 0- 22% baume 0- 25%Vol 0-40% brix optical refractometer	0-20%Baume 0- 25%Vol	0.2%Baume 0.2%Vol	±0.2%Baume ±0.2%Vol	Yes	10-30° C	Pure aluminium
YR05935	ATC Oe 0- 140°Oe 0- 27KMW (Babo)	0-140 °Oe 0-27 KMW (Babo) 0- 32%Brix (Sacch.)	±1°Oe ±0.2KMW ±0.2%Brix	±1°Oe ±0.2KMW ±0.2%Brix	Yes	10-30° C	Pure aluminium
YR05936	ATC Oe O- 190°Oe O- 38KMW (Babo)	0-190 °Oe 0-38 KMW (Babo) 0- 44%Brix (Sacch.)	±2°0e ±0.4KMW ±0.4%Brix	±2°0e ±0.4KMW ±0.4%Brix	Yes	10-30° C	Pure aluminium



#Letsgivemore Thanks to your purchase, a donat

All rights reserved © KALSTEIN France S. A. S.,

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S



www.kalstein.eu





#LETSGIVEMORE

WITH THE ADCQUISITION OF A KALSTEIN EQUIPMENT

YOU MAKE A CONTRIBUTION TO:





https://onetreeplanted.org





ANY DOUBT? Contact us!



Articles Kalstein (Social Networks and Official Websites)





PARIS - FRANCE **HEADQUARTER** 2 Rue Jean Lantier 75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89 **Cel:** +33 (0) 6 80 76 07 10 sales@kalstein.eu


https://kalstein.eu/ https://www.kalstein.fr/

Kalstein Training Videos (Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes 445 views • 1 week ago

