




# FURNITURE SECTOR

## ALL STAINLESS STEEL WORK BENCH



**Kalstein**®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



# ALL STAINLESS STEEL WORK BENCH

## KALSTEIN® All Stainless Steel Work Bench All Steel Work Bench Model YR06018 // YR06024

CE



#Let'sgivemore 

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S



**KALSTEIN® All Stainless Steel Work Bench**

**All Steel Work Bench**

**Model YR06018 // YR06024**



**Product Description:**

- ✓ Worktop: Chemical resistant phenolic resin, thickness 12.7 mm with edge double layer thickened to 25mm, Resistant to strong acid, strong alkali, anti-corrosion.
- ✓ Frame: Steel frame, epoxy resin powder spray on the surface, high temperature curing, high adhesion, featured with anti-scratch, anti-corrosion.
- ✓ Cabinet: 1.0mm thick cold rolled steel sheet coated with epoxy resin.
- ✓ Sideway: 3 section Sideway, no noise.
- ✓ Adjustable feet: Corrosion-resist high strength.
- ✓ Hinge: Open type high quality hinge.
- ✓ Handle: Stainless steel bright type handle.
- ✓ Standard color: Black worktop+ grayish cabinet.
- ✓ Socket: Not supply.



**Technical Specifications:**

Model	Description	Size mm (W*D*H)
YR06018	Side workbench w/ Reagent Rack (optional)	L*750*800 mm
YR06019	Central bench with double sides (one meter contain one seat area/leg room) w/ Reagent Rack (optional)	L*1500*800 mm
YR06020	Corner workbench	1000*1000*800 mm
YR06021	Sink workbench	950*750*800 mm
YR06022	Faucet and pp sink – high density PP	380*280*280 mm
YR06023	Drip rack- high density PP Type: single side	550*400*120 mm
YR06024	Balance bench (scale)	900*600*800
	Frame: Using 30*50*2 mm square steel, in order to strengthen the overall stability, square steel can be filled with sand. Epoxy resin powder sprayed for steel frame.	
	Worktop: 12.7 mm thickness Chemical resistant phenolic resin, put the balance position with 18 mm thick marble countertop, shock-reduction design, to ensure the accuracy use of laboratory equipment.	
	Adjustable caster: Use shock absorption casters of our own design.	
	Two sockets with anti-splash cover (two & three plugs)	
	Built-in three-stage shock absorption spring under main table: Highest accuracy millionth	



**#Letsgivemore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
 2 Rue Jean Lantier • 75001 Paris •  
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S. A. S**





 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación  
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale


<https://www.humatem.org/>



**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed

**The 3 articles of the week!**

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**1**  
Read this ARTICLE NOW



**How does the climate affect your microscope?**

**Kalstein Blog**

**Article of the week:**  
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**2**  
Read this ARTICLE NOW



**What is the importance of the laboratory refrigerator?**

**Kalstein Blog**

**Article of the week:**  
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**

2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89

**Cel:** +33 (0) 6 80 76 07 10

[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

