


# Kalstein A different accompaniment, at your service

# LABORATORY LINES BALANCES





All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu





OUR SERVICES

#### **Benefits and Support**

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



#### **Quick Response**

Our work team is always available to response all your consults or questions, in order to support you in any situation.





#### #Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



#### **Technical Support**

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.





#### **Delivery Logistics**

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or.air.



#### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



#### **Analytical Balances**

#### **Model Series YR05573**

#### Function:

- ✓ Internal auto calibration
- ✓ 5 Inch touch panel
- ✓ One key/Time trigger internal calibration
- ✓ Electromagnetic force sensor
- ✓ Aluminium shell/Large space wind proof cover
- ✓ Tare function/Counting/Density/Percentage
- ✓ Unit conversion(g/mg/ct/oz)...
- ✓ Under weighing
- ✓ Sensitivity/Speed set
- ✓ Overload alarm/Fault alarm/Level indicator
- ✓ Option: Density kit (Solid/Liquid)




































#### **Product application**











#### #Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu





### **Analytical Balances**

#### **Model Series YR05573**



Model	YR05573	YR05574	YR05575
Capacity	100g	210g	220g
Readability	0.0001g		
Min Weighing	0.0004g		
Repeatability	±0.0002g		
Linearity	±0.0003g		
Stable Time	<b>≤</b> 3\$		
Pan Size	Ø80mm		
Appearance Size	340x215x350mm		
<b>Draft Shield Size</b>	240x190x265mm		
Communication	RS232(Standard)/ USB(Optional)		
Temperature	15-35° <b>C</b>		
Calibration	Internal Calibration		
Cal.Weight	100g		
Option	Printer		
Packing Size	Inner box:415x285x430mm Outer box:475x315x485mm		
Gross Weight	11kg		



#### #Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu

KALSTEIN FRANCE, S. A. S







#LETSGIVEMORE >

# WITH THE ADCQUISITION OF A KALSTEIN EQUIPMENT

YOU MAKE A CONTRIBUTION TO:





https://onetreeplanted.org



Solidarité et coopération biomédicale

https://www.humatem.org/



https://www.maniapurefoundation.org/



## ANY DOUBT? Contact us!

#### **Articles Kalstein**

(Social Networks and Official Websites)







Article of the week: How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites





Article of the week: What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites https://kalstein.eu/articles



PARIS - FRANCE HEADQUARTER

2 Rue Jean Lantier 75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89 **Cel:** +33 (0) 6 80 76 07 10 sales@kalstein.eu

https://kalstein.eu/ https://www.kalstein.fr/



(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago

