




LABORATORY SECTOR

CLIMATIC CHAMBER



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Climatic Chamber
Artificial Climate Incubator
Model YR05328 // YR05330



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Climatic Chamber

Artificial Climate Incubator

Model YR05328 // YR05330



Bacteria Incubator:

Quality compressor with delay protection function; without fluorine refrigerant high efficient, more energy saving; built-in tempered glass viewing window and high -capacity intelligent humidification systems; unique air circulation to ensure uniform temperature. Simulate natural climate and environment, day and night, control temperature and humidity is an ideal environment test equipment for such as biology, medicine, agriculture and animal husbandry, forestry, environmental science research departments.

Suitable for plant cultivation and germination, seedling cultivation, bacteria, culture and preservation of microorganisms, small animals and insects feeding and so on.

Product Introduction:

Applies to plant seed germination, seedling cultivation, the cultivation and preservation of bacteria, microbes, small animals, insects, feeding, product quality testing and other uses of light, constant temperature experiments. Simulate natural climate and environment, through the light control simulation during the day and night, the humidity inside (this feature is limited to artificial climate incubator, the temperature can be adjusted to a constant temperature degree is an ideal environment test equipment for such as biology, medicine, agriculture and animal husbandry, forestry, environmental science research departments.

Features:

- ✓ Balanced refrigeration, small temperature fluctuations, and set the automatic defrost, multiple layers of security protection, required to meet the long-running, with the well-known international brand compressor with delay protection function and smooth running.
- ✓ Introduction of green philosophy, the use of Freon R134 refrigerant and make the system more energy efficient, healthier environment.
- ✓ The temperature control system with upgraded LCD display, programmable microcomputer intelligent PID control technology, with timing, ultra-high temperature alarm function.
- ✓ The unique air circulation, soft wind, to ensure temperature uniformity.
- ✓ Illumination five level can be adjusted, liner is made of mirror stainless steel, up and down pitch adjustable shelf.
- ✓ The double door design, built-in tempered glass viewing window, the items inside can be clearly observed.
- ✓ Artificial climate chamber with built-in high-capacity intelligent humidification systems, and automatic replenishment by the external tank.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Climatic Chamber Artificial Climate Incubator Model YR05328 // YR05330



Technical Parameters:

Product Name	Illumination (Ray radiation) incubator YR series		
Model	YR05328	YR05329	YR05330
Temp range °C	With illumination 10~60°C, With no illumination 5~60°C		
Humidity range	50%~90% humidity fluctuation±5%RH		
Temperature uniformity values °C	±1		±1.5
Illumination LX	0~10000	0~15000	0~15000
Illumination set range	Five level (20%、40%、60%、80%、100%)		
Timing range	0~99h		
Power W	900	1000	1200
Voltage V	220		
Accuracy °C	±0.1		
Volume L	150	250	350
Inner size (mm) (H*W*D)	800*500*380	1000*550*450	1000*580*600
Overall dimension	1525*735*680	1775*785*750	1775*815*900
(mm)			
(H*W*D)			
Packing size (mm)	1680*810*760	1925*860*830	1925*890*980
(H*W*D)			
Net weight (KG)	120	155	188
Shipping weight (KG)	128	165	198



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:

 **Fundación
Jacinto Convit**
<https://www.jacintoconvit.org/>


 **ONE
TREE
PLANTED**

<https://onetreepanted.org>

 **Humatem**
Solidarité et coopération
biomédicale
<https://www.humatem.org/>

 **MANIAPURE
FOUNDATION USA**
<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

1
Read this ARTICLE NOW



How does the climate affect your microscope?

microscope is an instrument that allows observing and measuring small objects. This is achieved through a series of lenses that focus light on the object being observed.

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

2
Read this ARTICLE NOW



What is the importance of the laboratory refrigerator?

Refrigerator is a device that is one of the most important equipment that helps to maintain, in a controlled environment (refrigerated space), various fluids and substances.

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

