




# LABORATORY LINES

## pH METERS



# Kalstein®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.





## KALSTEIN® Benchtop pH Meter Model YR01824



### YR01824 Benchtop pH Meter

pH electrode, pH buffer pouches [pH4.01, 7.00, 10.01], power adapters.

#### Features:

- ✓ 2 point calibration with automatic recognition for USA and NIST buffers
- ✓ Automatic electrode diagnosis shows the pH slope
- ✓ Manual temperature compensation ensures accurate readings over the entire range

#### Applications:

- ✓ Education: Chemical experiments
- ✓ Laboratory: Routine measurements and quality controls

#### Ordering Information:

- ✓ YR01824 pH meter
- ✓ E201-BNC pH electrode
- ✓ pH4.01/7.00/10.01 buffer reagents Electrode holder
- ✓ DC9V power adapter



**#Let'sgive more** ❤️  
Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
Optimum Business Center 450 Rue Baden Powell,  
• 34000 Montpellier, France.  
Tit: +33 467158849 / +33 680760710 / +33 663810023  
<https://kalstein.eu>  
KALSTEIN FRANCE, S.A. S



## KALSTEIN® Benchtop pH Meter

Model YR01824

	Model	YR01824
pH	Range	0.00-14.00pH
	Accuracy	±0.05pH
	Resolution	0.01pH
	Calibration Points	2
	pH Buffer Options	USA (pH4.01/7.00/10.01) or NIST (pH4.01/6.86/9.18)
mV	Range	±1999mV
	Accuracy	±10mV
	Resolution	1mV
Other Functions	Temperature Compensation	0~100°C, manual
	Connector	BNC
	Display	5.5" custom LCD
	Power	DC9V, using AC adapter, 220V/50Hz
	Dimensions	210(L)×205(W)×75(H)mm
	Weight	1.5KG



#Letsgivemore

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
Optimum Business Center 450 Rue Baden Powell,  
• 34000 Montpellier, France.  
Tél: +33 467158849 / +33 680760710/ +33 663810023  
<https://kalstein.eu>  
KALSTEIN FRANCE, S.A. S







 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación  
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale

<https://www.humatem.org/>




**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>



**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed  
**The 3 articles of the week!**  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW  
**How does the climate affect your microscope?**  
A microscope is an instrument that allows observing and measuring objects that are too small to be seen with the naked eye. This is achieved through a complex system of lenses and mirrors, which form...

**Kalstein Blog**  
**Article of the week:**  
How does the climate affect your microscope?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW  
**What is the importance of the laboratory refrigerator?**  
A refrigerator in a laboratory is one of the most important pieces of equipment. It is used to maintain, in a controlled environment (refrigerated space), various fluids and substances...

**Kalstein Blog**  
**Article of the week:**  
What is the importance of the laboratory refrigerator?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



MONTPELLIER - FRANCE  
**HEADQUARTER**  
Optimum Business Center,  
450 Rue Baden Powell,  
34000 Montpellier - France.

**Fax:** +33 (0) 4 67 15 88 49  
**Tif:** +33 (0) 6 80 76 07 10  
sales@kalstein.eu

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

**Official Websites for Latin America.**

