



LABORATORY LINES

MICROSCOPES



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Microscopes Binocular Student Biological Microscope

Model YR0242 // YR0242-1



Product Description

1. Anti-mould System Head Available
2. 195mm Achromatic Objectives, Mechanical Tube Length 160mm
3. WF10x/18mm, WF16x/11mm Wide Field Eyepiece
4. Monocular, Binocular, Trinocular Head Optional
5. Low Position Coaxial Coarse & Fine Focus System
6. Double Layer Mechanical Stage



Model	YR0242	YR0242-1
Magnification	40x~1000x	
Head	Sliding Binocular Head, 45° Inclined	0
	Sliding Binocular Head, 45° Inclined, Anti-Mould System	0
Eyepiece	Wide Field WF10x/18mm	
Nosepiece	Quadruple Nosepiece	
Objective	Achromatic 4x/0.10, 10x/0.25, 40xs/0.65, 100xs/1.25(Oil)	
Stage	Double Layer Mechanical Stage, 140*140mm, Mvoing Range 50*75mm	
Focusing	Coaxial Coarse & Fine Focusing. Fine Focusing Scale 0.002mm, Coarse Focusing Stroke 37.7mm/Rotation, With Tension Adjustment, With Up Limit	
	Rubber Cover On Focusing Knob	0
Condenser	N.A.1.25 Abbe Condenser, Iris Diaphargm, With Filter Holder	
Filter	Blue, Yellow, Green	
Light Source	Plan-Concave Mirror	
	Built-in Illumination 6V20W Halogen Lamp, 100~240V	
YR0242 Biological Microscope Optional Accessories		
Head	Sliding Trinocular Head, 45° Inclined, Anti-Mould System	YR0242-1-op1
	Seidentopf Binocular Head, 30° Inclined, Anti-Mould System	YR0242-1-op2
	Seidentopf Trinocular Head, 30° Inclined, Anti-Mould System	YR0242-1-op3
Eyepiece	Wide Field Eyepiece P16x/11mm	YR0242-1-op4
Nosepiece	Quintuple Nosepiece	YR0242-1-op5



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
A microscope is an instrument that allows observing and measuring objects at a very small scale. This is achieved through a complex system of lenses and mirrors, which form...

Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
A refrigerator is a device that is one of the most common pieces of equipment used in laboratories. It is a controlled environment (refrigerated space), various fluids and substances...

Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

