

# **LABORATORY SECTOR**

REFRACTOMETER







All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S





# OUR SERVICES Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



#### Online Inductions and Trainings In any part of the world, receive your

induction or training from our specialized team of engineers



#### **Quick Response** Our work team is always available to response all your consults or questions, in order to support you in any situation.





## #Letsgivemore 💙

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



#### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



## **Delivery Logistics**

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or.air.



#### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.





## KALSTEIN® Refractometer Brix Refractometers Model YR05872 // YR05884





#Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.

All rights reserved © KALSTEIN France S. A. S.,



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S



www.kalstein.eu

#### KALSTEIN® Refractometer Brix Refractometers Model YR05872 // YR05884



## Product Description:

Brix Refractometer is frequently used by grape growers, vegetables growers, home brewers as well as food producers for Brix testing of sucrose levels in their fruit. Brix refers to the measurement of the percentage of sugars and dissolved solids in a given weight of fresh plant juice. These solids include sucrose, fructose, vitamins, minerals, amino acids, proteins, hormones, and other solids commonly found in a plant's juice. This RHB-05 is specially designed for cutting liquids testing. Also, it can be used to test the low brix solution.

## Features:

- ✓ High quality body construction
- Uses ambient light only which means battery or power source is not required
- $\checkmark$  Easy to focus and calibrate
- Cushioned with soft & comfortable non-slip rubber
- ✓ With ATC function(Auto Temperature Compensation)
- ✓ Durable and built to last long

All rights reserved © KALSTEIN France S. A. S.,





**#Letsgivemore** Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S. 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S



## KALSTEIN® Refractometer Brix Refractometers Model YR05872 // YR05884



#### **Technical Specifications** :

| Model   | Description  | Measurement<br>Range | Min. Div. | Accuracy    | АТС | ATC temp. | Material      |
|---------|--|----------------------|-----------|-------------|-----|-----------|---------------|
| YR05872 | ATC Brix 0-5%<br>optical<br>refractometer                                  | 0-5%Brix             | 0.1%Brix  | ± 0.10%Brix | Yes | 10-30° C  | Pure aluminum |
| YR05873 | ATC Brix 0-10%<br>optical<br>refractometer                                 | 0-10%Brix            | 0.1%Brix  | ± 0.10%Brix | Yes | 10-30° C  | Pure aluminum |
| YR05874 | ATC Brix 0-18%<br>optical<br>refractometer                                 | 0-18%Brix            | 0.1%Brix  | ± 0.10%Brix | Yes | 10-30° C  | Pure aluminum |
| YR05875 | ATC Brix 0-20%<br>optical<br>refractometer                                 | 0-20%Brix            | 0.1%Brix  | ± 0.10%Brix | Yes | 10-30° C  | Pure aluminum |
| YR05876 | ATC Brix 0-32%<br>optical<br>refractometer                                 | 0-32%Brix            | 0.2%Brix  | ± 0.20%Brix | Yes | 10-30° C  | Pure aluminum |
| YR05877 | ATC Brix 0-50%<br>optical<br>refractometer                                 | 0-50%Brix            | 0.5%Brix  | ± 0.50%Brix | Yes | 10-30° C  | Pure aluminum |
| YR05878 | ATC Brix 0-55%<br>optical<br>refractometer                                 | 0-55%Brix            | 0.5%Brix  | ± 0.50%Brix | Yes | 10-30° C  | Pure aluminum |
| YR05879 | ATC Brix 28-<br>62% optical<br>refractometer                               | 28-62%Brix           | 0.2%Brix  | ± 0.20%Brix | Yes | 10-30° C  | Pure aluminum |
| YR05880 | ATC Brix 45-<br>82% optical<br>refractometer                               | 45-82%Brix           | 0.5%Brix  | ± 0.50%Brix | Yes | 10-30° C  | Pure aluminum |
| YR05881 | ATC Brix 58-<br>92% optical<br>refractometer                               | 58-92%Brix           | 0.5%Brix  | ± 0.50%Brix | Yes | 10-30° C  | Pure aluminum |
| YR05882 | ATC Brix 0-80%<br>optical<br>refractometer                                 | 0-80%Brix            | 1%Brix    | ± 1.0%Brix  | Yes | 10-30° C  | Pure aluminum |
| YR05883 | ATC Brix 0-90%<br>optical<br>refractometer                                 | 0-90%Brix            | 1%Brix    | ± 1.0%Brix  | Yes | 10-30° C  | Pure aluminum |
| YR05884 | ATC 58-90%Brix<br>38-43Be' 12-<br>27%water LED<br>optical<br>Refractometer | 0-90%Brix            | 1%Brix    | ± 1.0%Brix  | Yes | 10-30° C  | Pure aluminum |



**#Letsgivemore** Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S



#### www.kalstein.eu





#LETSGIVEMORE

## WITH THE ADCQUISITION OF A KALSTEIN EQUIPMENT

YOU MAKE A CONTRIBUTION TO:





https://onetreeplanted.org





## ANY DOUBT? Contact us!



#### Articles Kalstein (Social Networks and Official Websites)





PARIS - FRANCE **HEADQUARTER** 2 Rue Jean Lantier 75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89 **Cel:** +33 (0) 6 80 76 07 10 sales@kalstein.eu

https://kalstein.eu/ https://www.kalstein.fr/

#### Kalstein Training Videos (Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes 445 views • 1 week ago

