




MEDICAL LINES

PATIENT MONITOR



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.





YR02190A Central Monitoring System

YR02190A Central Monitoring System provides the LAN or wireless LAN networking solution for the hospital departments. Combining the flexible monitoring with useful clinical tools and intuitive user interface layouts, **YR02190A** cooperates perfectly with all **Kalstein France** patient monitors. It improves the efficiency by helping you care more on patients.



#Letsgivemore ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tit: +33 467158849 / +33 680760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



Feature

- Enable to monitor 32 beds and be expandable to 64 beds in one system
- Include clinical support tools: real-time trend display, tabular and graphic trend review, etc.
- Improve the patient care efficiency by the remote control

Technical Specifications

Recommend Hardware Configuration

OS	Windows XP or Windows 7
CPU	Intel Core Dual, 3.0G Hz
RAM	2 GB
Display	19", 1280 x 1024 resolution
Graphic Memory	256 MB
Hard disk space	≥ 250GB
Internet connection speed	≥ 10M bps

Environment

Operating Temperature	5 ~40 °C
Storage Temperature	-20~55 °C
Relative Humidity	≤ 80 %
Atmosphere pressure	86~106 kPa
Power	AC 200 V/50 Hz

Safety

Meet GB9706.15-2008 Medical electrical equipment

Real-time View

2-channel ECG waveform
1-channel RESP waveform
1-channel waveform PLETH
HR, RR in ECG, Oxygen saturation,
PR, Systolic, diastolic and mean
pressure in NIBP, 2-Temp, ST value,
EtCO2, InCO2, AWRR, 2-IBP

Review

NIBP review for each monitor
Trend review for each monitor
Alarm review for each monitor
History patient monitoring information review
Full disclosure ECG waveform review
Review information output (by printer)



#Let'sgive more ❤️


Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Titl: +33 467158849 / +33 680760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
A microscope is an instrument that allows observing and measuring objects that are too small to be seen with the naked eye. This is achieved through a complex system of lenses and mirrors, which form an image of the object being observed.
Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
A refrigerator in a laboratory is one of the most important pieces of equipment. It is used to maintain, in a controlled environment (refrigerated space), various fluids and substances.
Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



MONTPELLIER - FRANCE
HEADQUARTER
Optimum Business Center,
450 Rue Baden Powell,
34000 Montpellier - France.

Fax: +33 (0) 4 67 15 88 49
Tif: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.net/>
<https://www.kalstein.fr/>

Official Websites for Latin America.

