




LABORATORY LINES

COAGULATION ANALYZER



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.




KALSTEIN® Coagulation Analyzer Model YR05117



Features

- ✓ Double testing channels, 2 different items analysis in the same time
- ✓ Special testing cup location system, sample bracket
- ✓ Incubate precision 37°C±0.3°C
- ✓ Up to 10.000 testing result restored ,Inner thermal-sensitive printer
- ✓ Electronically-linked pipet optional
- ✓ Perform all routine tests such as PT, APTT, TT, FiB and clotting factors
- ✓ RS232 interface, PC connecting



#Let'sgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



KALSTEIN® Coagulation Analyzer Model YR05117



Test channel number	2 channel;2 different parameters analysis in the same time
Sample incubator	16positions
Reagent incubator block	4 positions
Software timer	2
Test wavelength	470nm
Memory	10000 test result
Temp control precision	37±0.3°C
Sample consume	20uL-40uL
Reagent consume	20uL-40uL
Screen	240*128 LCD
Coincidence indicator	≤3%
Notation of Result	S.%,PTR,INR,g/L
Operation Language	English
Power supply	220V±15%, 50-60Hz
Dimension	395mm (L)*320mm (W)*135mm(H)
Weight	6 KG



#Letsgivingore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tit: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



KALSTEIN® Coagulation Analyzer Model YR05118



Features

- ✓ Low reagent consumption less than 20ul, open reagent.
- ✓ Special testing cup location system, adding sample bracket
- ✓ Advanced optical light matching system ensure precision result.
- ✓ Optional PT Derived method for testing fibrinogen
- ✓ Decrease the dosage of reagent
- ✓ Top quality thermostat, precision 37°C ± 0.3 °C
- ✓ Dispense precision: CV <2%
- ✓ Memory for 10.000 testing result
- ✓ Internal thermal-sensitive printer
- ✓ Electronically-linked pipette optional
- ✓ Perform all routine tests such as PT.APTT.TT Fibrinogen and clotting factors
- ✓ RS232 Interface, PC Connecting

Usual text list:

- ✓ Prothrombin PT
- ✓ Activated Partial thromboplastin time APTT
- ✓ Thrombin Time TT
- ✓ Fibrinogen FIB
- ✓ RT;VT;VIII;IX;XI;XII;II;VII;X;PS;PC;Xa;etc



#Letsgivemore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



KALSTEIN® Coagulation Analyzer Model YR05118



Model	YR05118
Test Channel Number	4 Channel; 4 Different Parameters analysis in the same time
Sample Incubation Block	24 Positions
Reagent Incubation Block	6 Positions
Software Timer	4
Test Wavelength	470 nm
Memory	10 000 test result
Temperature Control Precision	37 ± 0.3°C
Sample Consume	20 uL – 40 uL
Reagent Consume	20 uL – 40 uL
Coincidence Indicator	≤3%
Notation of Result	S, %, PTR, INR, g/L
Operation Language	English
Screen	240 x 128 LCD
Power Supply	220 VAC ± 15% 50-60Hz
Dimensions (LxWxH)	395 x 320 x 135 mm
Work Environment	15 °C – 30°C wet ≤90%
Weight	7 Kg



#Letsgivingmore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:

 **Fundación
Jacinto Convit**
<https://www.jacintoconvit.org/>


 **ONE
TREE
PLANTED**

<https://onetreepanted.org>

 **Humatem**
Solidarité et coopération
biomédicale
<https://www.humatem.org/>

 **MANIAPURE
FOUNDATION USA**
<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
microscope is an instrument that allows observing and measuring small objects. This is possible through a series of lenses that focus light from the object being observed.

Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
Refrigerator is a device that is one of the most common equipment that is used to maintain, in a controlled environment (refrigerated space), various fluids and substances.

Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



MONTPELLIER - FRANCE
HEADQUARTER
Optimum Business Center,
450 Rue Baden Powell,
34000 Montpellier - France.

Fax: +33 (0) 4 67 15 88 49
Tel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

Official Websites for Latin America.

