




# LABORATORY LINES

COLOR LIGHT BOX



## Kalstein®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



## KALSTEIN® Color Light Box Color Viewer Light Box with Drawer Model YR05545//YR05547



### Product Description

Model YR05545 // YR05547 color viewer light table features with upper and lower light and three drawers. It has 3 drawers more than YR05545 // YR05547, and others are the same with YR05545 // YR05547. It can choose one or two or three light sources from D65, D50, TL84 in the upper part, which will have 4 pcs installation there. And it also can choose only one light source from D65, D50 or TL84 in the lower part, which is with 6 pcs installation there. Different light source will have different price.

We offer standard light sources commonly used in the printing industry to meet ISO 3664:2009 standards.

YR05545 // YR05547 vertical standard light source color viewer is suitable for the printing workshop or proofing center with high environmental requirements. It can be applied to the printing factory with beautiful appearance, offset printing, flexo printing and gravure printing. It is also directly able to see the reflective manuscript.

YR05545 // YR05547 series vertical sample table is equipped with frequently-used light sources in printing industry, can meet the needs of ISO 3684:2009. It is also widely used for color checking in the field of printing and dyeing, textiles, plastics, hardware etc, which will provide standard environment for quality checking in the process of sample comparison, production, quality inspection, acceptance etc. That is to say, it can help user improve the product quality and market competition.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
Optimum Business Center 450 Rue Baden Powell,  
• 34000 Montpellier, France.  
Tél: +33 467158849 / +33 680760710/ +33 663810023  
<https://kalstein.eu>  
KALSTEIN FRANCE, S.A. S



## KALSTEIN® Color Light Box Color Viewer Light Box with Drawer

Model YR05545//YR05547



### PRODUCT FEATURES

Light Sources: can meet the international standard of ISO 3664

Lamps type: equipped with special Graphic lamps for import printing machine

Color-rendering index: CRI=98

Worktable material: metal plate desk, firm & durable

Operation desk: is incorporated with an intake ramp designed to meet the human visual observation

Hanging light sources: aluminium alloy stuff and streamline from design, fashion, style. Import reflector, uniform light cover will ensure the illuminating uniformity

### APPLICATION INDUSTRY

YR series vertical sample table is equipped with frequently-used light sources in printing industry, can meet the needs of ISO 3664:2009. It is also widely used for color checking in the field of printing and dyeing, textiles, plastics, hardware etc, which will provide standard environment for quality checking in the process of sample comparison, production, quality inspection, acceptance etc. That is to say, it can help user improve the product quality and market competition.



Printing Shop



Printing & Dyeing



Textiles



Plastics



Hardware



Others



#Letsgivemore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
Optimum Business Center 450 Rue Baden Powell,  
• 34000 Montpellier, France.  
Tél: +33 467158849 / +33 680760710/ +33 663810023  
<https://kalstein.eu>  
KALSTEIN FRANCE, S.A. S



## KALSTEIN® Color Light Box Color Viewer Light Box with Drawer

Model YR05545//YR05547



### PRODUCT PARAMETERS

| Upper light source with three drawers       |                                |   |                 |                   |            |             |
|---|--------------------------------|---|-----------------|-------------------|------------|-------------|
| The upper configured with 3 groups of light | Model and specification        | standard light source                               |                 | Color Temperature | Power      | Quantity    |
|   | Vertical viewing table YR05545 | D65 single light source, with three drawers         |                 | 6500K             | 36W        | 1 x 4 Tubes |
|   | Vertical viewing table YR05546 | D65 D50 dual light source, with three drawers       |                 | 6500 5000K        | 36W        | 2 x 4 Tubes |
|   | Vertical viewing table YR05547 | D65 D50 TL84 Three light sources,with three drawers |                 | 6500 5000 4000K   | 36W        | 3 x 4 Tubes |
| Lower light source with three drawers       |                                |   |                 |                   |            |             |
| Optional standard light source              | Vertical viewing table YR05545 | D65 D50 TL84 single light source                    | 6500 5000 4000K | 36W               | 1 x 6Tubes |             |
| Basic Parameters                            |                                |   |                 |                   |            |             |
| Power Supply                                |                                | AC220V 50HZ   |                 |                   |            |             |
| Workable size                               |                                | 1350*900(L*W)                                       |                 |                   |            |             |
| Vertical size                               |                                | 1350*2174(L*W)                                      |                 |                   |            |             |
| External size                               |                                | 1350*900*217(L*W*H)                                 |                 |                   |            |             |
| Drawer Height(mm)                           |                                | 120*3 pcs   |                 |                   |            |             |



#Letsgivemore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
Optimum Business Center 450 Rue Baden Powell,  
• 34000 Montpellier, France.  
Tél: +33 467158849 / +33 680760710/ +33 663810023  
<https://kalstein.eu>  
KALSTEIN FRANCE, S.A. S







#LETSGIVEMORE

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:




<https://onetreepanted.org>





**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed

**The 3 articles of the week!**

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW

**How does the climate affect your microscope?**

**Kalstein Blog**

**Article of the week:**  
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW

**What is the importance of the laboratory refrigerator?**

**Kalstein Blog**

**Article of the week:**  
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**

2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89  
**Cel:** +33 (0) 6 80 76 07 10  
[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

