



LABORATORY LINES

VISCOMETERS



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Viscometers

Digital Rotating Viscometer Brookfield Viscosity Tester

Model YR05861 // YR05863



Product Description

1. Test range from 1 to 2,000,000 mpa.s (cps)
2. High accuracy
3. Automatic detection
4. Easy to operate



Technical Specifications

Model	YR05861	YR05862	YR05863
Measurement Range	1-2,000,000 mPa·S		
Display	LCD		
Spindles Specification	0, 1, 2, 3, 4	0, 1, 2, 3, 4, 21, 27, 28, 29	
Rotating Speed of Spindles	0.3/0.6/1.5/3/6/12/30/60	0.3/0.6/1.5/3/6/12/30/60	
	0.5/1/2/2.5/4/5/10/20/50/100	0.5/1/2/2.5/4/5/10/20/50/100	
	0	0.5/1/2/5/10/20/50	
Interface	RS232		
Automatic Switch	free selection of proper rotating speed or Spindles number.		
Stable Cursor	Strip-type cursor fully charged displayed when reading basically stable		
Measurement accuracy	±1%		
Measurement Repeatability	±0.3%		
Viscosity Resolution	1 mPa·S (If with No. 0 spindles, the resolution is 0.1 mPa.s)		
Temperature probe	-50 to 99 degree		
Temperature precision	± 0.5 degree		



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



KALSTEIN® Viscometers

Digital Rotating Viscometer Brookfield Viscosity Tester

Model YR05861 // YR05863



Technical Specifications

Model	YR05861	YR05862	YR05863
Temperature resolution	0.1 degree		
Power	AC 220V ± 10% 50 Hz ± 10%		
Working Condition	Temperature 5-35 degree / Relative humidity no more than 80%		
Attachments	No. 0, 1, 2, 3, 4 Spindles	No. 0, 1, 2, 3, 4, 21, 27, 28, 29 Spindles	
		Temperature probe	
		Software	
		Adapter	
		Foundation Plate	
		Up-Down Frame	
Additional Information	English Manual		
	0.3/0.6/1.5/3/6/12/30/60 rpm/min	0.3/0.6/1.5/3/6/12/30/60 rpm/min	0.5/1/2/5/10/20/50 rpm/min
Mpa.S	1-2,000,000 cps		



#Letsgivingore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tit: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



KALSTEIN® Viscometers

Digital Rotating Viscometer Brookfield Viscosity Tester

Model YR05861 // YR05863



Features and accessories



Spindles



Packing: 42*34*35cm/carton



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
A microscope is an instrument that allows observing and measuring objects at a very small scale. This is achieved through a complex system of lenses and mirrors, which form...

Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>




2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
A refrigerator is a device that is one of the most common pieces of equipment used in laboratories. It is a controlled environment (refrigerated space), various fluids and substances...

Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

