




# LABORATORY LINES

## DRY BATH INCUBATOR



**Kalstein**®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



## Dry Bath Incubator



YR0A150-1



YR0A150-1A



YR0A150-2



YR0A150-4

### Introduction

**Kalstein** digital dry bath is controlled by microcomputer, widely used in preservation and reaction of the samples, DNA amplification and pre-denaturation of electrophoresis, serum coagulation, etc.

### Features

- LED display with timer
- Automatic fault detection with buzzer alarm
- Temperature calibration function
- Various blocks, easy to clean and disinfect
- Built-in over-high temperature protection



Various blocks



High precision



Over-temperature protect



Timer

### Product parameters

Model	YR0A150-1	YR0A150-1A	YR0A150-2	YR0A150-4	Model	YR0A150-1	YR0A150-1A	YR0A150-2	YR0A150-4
Temperature Control Range	RT+5°C ~150°C				Display accuracy	0.1°C			
Heating Up Rate	≤ 10M (20°C to 100°C)				Timer	0~99h59m			
	≤ 30M (20°C to 150°C)								
Temperature Stability @40°C	±0.3°C				Max.Temperature	150°C			
Temperature Stability @100°C	±0.5°C				Power	200W	400W	600W	
Temperature Stability @120°C	±1°C				Dimension (W×D×H)	220*260*100 (mm)			220*360*100 (mm)



### #Let's Give More

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S





“D1




“D2



“D3



“D4



“D5



“D6



“D7



“D8



“D9



“D10



“D11



“D12



“D13



“D14



“D15



“D16



“D17



“D18

Important Notes: D17/D18 are used in YR0A1501A-1A

“Optional Blocks							
Type	Tube diameter	Capacity	size (W×D×H)	Type	Tube diameter	Capacity	size (W×D×H)
D1	6mm	48	96.5×76.5×50 (mm)	D10	26mm	8	96.5×76.5×50 (mm)
D2	7mm	48	96.5×76.5×50 (mm)	D11	28mm	6	96.5×76.5×50 (mm)
D3	10mm	24	96.5×76.5×50 (mm)	D12	38mm	3	96.5×76.5×50 (mm)
D4	12mm	24	96.5×76.5×50 (mm)	D13	0.2ml tubes	48	96.5×76.5×50 (mm)
D5	13mm	24	96.5×76.5×50 (mm)	D14	0.5ml tubes	48	96.5×76.5×50 (mm)
D6	15mm	12	96.5×76.5×50 (mm)	D15	1.5ml tubes	24	96.5×76.5×50 (mm)
D7	16mm	12	96.5×76.5×50 (mm)	D16	2ml tubes	24	96.5×76.5×50 (mm)
D8	19mm	12	96.5×76.5×50 (mm)	D17	0.2ml	96	96.5×153×33.5 (mm)
D9	20mm	12	96.5×76.5×50 (mm)	D18	Flat	96 Elisa	96.5×153×22.5 (mm)



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
 2 Rue Jean Lantier • 75001 Paris •  
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
 KALSTEIN FRANCE, S. A. S







 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación  
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale

<https://www.humatem.org/>




**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>



**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed

**The 3 articles of the week!**

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW

**How does the climate affect your microscope?**

microscope is an instrument that allows observing and measuring objects that are too small to be seen with the naked eye. This is achieved through a complex system of lenses and mirrors, which form

**Kalstein Blog**

**Article of the week:**  
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW

**What is the importance of the laboratory refrigerator?**

refrigerator is a device that is one of the most common pieces of laboratory equipment. Its purpose is to maintain, in a controlled environment (refrigerated space), various fluids and substances

**Kalstein Blog**

**Article of the week:**  
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



**Kalstein France - Jacinto Convit Foundation / Alliances**  
23 views • 3 years ago



**Kalstein's Microscopes**  
445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**  
2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89  
**Cel:** +33 (0) 6 80 76 07 10  
[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

