




DENTAL SECTOR

DENTAL UNIT



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Dental Units
Efficient Work Dental Unit
Model YR06009



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Dental Units

Efficient Word Dental Unit

Model YR06009



Product Description:

Dental comprehensive treatment chair, Perfect chair back design, healthy and comfortable treatment space, bring upgrade experience to patients. New spring super sensitive touch screen, multi-function, fast and efficient operation. The YR06009 dental comprehensive treatment chair with interlocking function of hand piece and chair and anti-collision function of chair back, can prevent accidents in the treatment process and provide a sense of safety for doctors and patients. Perfect chair back design, healthy and comfortable treatment space, bring upgrade experience to patients. New spring super sensitive touch screen, multi-function, fast and efficient operation.



Features:

- ✓ The dental comprehensive treatment chair, New spring super sensitive touch screen, multi-function, fast and efficient operation. With interlocking function of headpiece and chair and anti-collision function of chair back, can prevent accidents in the treatment process and provide a sense of safety for doctors and patients. Perfect chair back design, healthy and comfortable treatment space, bring upgrade experience to patients.
- ✓ The lowest chair position can be as low as 380mm, which is convenient for all ages to see a doctor and ensures the best comfortable treatment posture.
- ✓ The doctor's chair is wide, ergonomic, and can be moved flexibly with a silent chair wheel, making the work easier.
- ✓ The 45 ° doctor and assistant operation panel is customer-oriented, with sensitive touch keys and convenient operation. It has three memory functions and other functions. Spitting can be reset with one button, which is simple, fast and efficient.
- ✓ Suction handle with regulator, compatible with central suction system. Optional top mounted model. The patient chair is made of super fiber leather, which is breathable, comfortable and has long service life.
- ✓ The precise four-way foot control and the one button flushing and water supply button liberate the doctor's hands. The back of the chair is linked with the patient's chair, and the ergonomic cushion makes the patient comfortable. Wide cast aluminum chassis , more stable and safer. Flexible swiveling and compact workflow assistant unit easy to reach, provide more space for assistant working closely with dentist and patient all the time.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Dental Units
Efficient Word Dental Unit
Model YR06009



Standard Configuration:

Model	YR06009
3-way syringe(cold/hot)	2sets
Powerful suction device	1set
Saliva ejector	1set
Dental light	1set
Warm water system	1set
Pure water supply system	1set
Rotatable ceramic spittoon	1set
DC 24V motor driven chair	1set
Multifunctional assisant device	1set
LED x-ray viewer	1set



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

