




# MEDICAL SECTOR

## WHEELCHAIR



**Kalstein**®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



**KALSTEIN® Wheelchair**  
**Electric Climb Stairs Wheelchair**  
**Model YR05446**



**#Letsgivemore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S



**KALSTEIN® Wheelchair**  
**Electric Climb Stairs Wheelchair**  
**Model YR05446**



**Product Description:**

Electric climb stairs wheelchair YR05446, turns dreams into reality; with this unique device, you can explore the world like never before and bring a new taste to your everyday routines! The YR05446 wheelchair can help you climb stairs, overcome obstacles, and travel in any direction by using the one joystick on the side of the wheelchair. They are extremely popular because they're easy to use and spare the users from having to use their arms to propel the chair. It also makes them less dependent on another person. These chairs are meant for everyone and can be especially useful for people with more serious disabilities who can't propel a manual chair. The newest designs come with some amazing features such as a speech device, various controls, options to adapt the body position, and more.

**Description:**

Electric climb stairs wheelchair YR05446, long-range, lithium wheelchair is the best stair-climbing wheelchair overall. It comes with a big lithium battery and is three times faster than other power wheelchairs. It also works well on all kinds of stairs regardless of the texture.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
 2 Rue Jean Lantier • 75001 Paris •  
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S. A. S**



## KALSTEIN® Wheelchair Electric Climb Stairs Wheelchair Model YR05446



### Technical Description:

<b>Model</b>	YR05446
<b>Frame</b>	Aluminum alloy
<b>Battery</b>	18Ah lithium battery
<b>Wheel size</b>	8 & 12 inches
<b>Net weight</b>	38kg
<b>Max loading</b>	150kg
<b>Climb/Descent</b>	three-speed
<b>Climbing ability</b>	500 stair steps
<b>Stair width required</b>	> 700mm
<b>Stair slope required</b>	<45 °
<b>Stair depth required</b>	1100mm
<b>Product size</b>	98*65*125cm
<b>Folding size</b>	110*65*35cm
<b>Seat width</b>	43cm
<b>Seat depth</b>	45cm
<b>Cushion</b>	Breathable sponge
<b>Color</b>	Black/Yellow/Red



#Letsgivemore ❤️

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S





 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación  
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale


<https://www.humatem.org/>



**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed

**The 3 articles of the week!**

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW

**How does the climate affect your microscope?**

**Kalstein Blog**

**Article of the week:**  
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW

**What is the importance of the laboratory refrigerator?**

**Kalstein Blog**

**Article of the week:**  
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**

2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89

**Cel:** +33 (0) 6 80 76 07 10  
[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

