




LABORATORY LINES

ELECTROPHORESIS



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



Mini-Protean Vertical Electrophoresis Cell

Model YR03425



Vertical gel electrophoresis is a more complex set-up compared to horizontal gel system. Researchers normally use this system to separate proteins instead of nucleic acids. However, before you separate proteins, you need to disrupt the quaternary structure of proteins into the linear strand

Characteristic features:

- ✓ The installation of the glass plate in just 15 seconds, and does not need the knob, fast and convenient.
- ✓ With a gel casting in original position, no need to move glasses from gel making to electrophoresis running. It is more convenient to observe gel preparing through both sides of the glasses.
- ✓ High transparent polycarbonate material injection moulding, impact resistance, high temperature resistance, corrosion resistance.
- ✓ Safety cover open button design, convenient opening the cover.
- ✓ With the background colour, easy to be observed in the sample adding and electrophoresis process.
- ✓ Abundant buffer not only guarantees cooling effects, but also keeps the pH Value stable during the whole experiment process.
- ✓ Auto-switch-off when the lid is removed.



#Letsgivemore

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tit: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



Mini-Protean Vertical Electrophoresis Cell

Model YR03425



Model	YR03425
Glass Plate Size (WxL)	100×100 mm
Gel Size (WxL)	82×88 mm
Gel Thickness	0.75 、 1.0、 1.5 (mm)
Gel Amount	1~2 (pieces)
Sample Volume	(0.75mm) 11、 15 wells; (1.0mm) 11、 15wells; (1.5mm) 11、 15wells
Buffer Volume	~750 (ml)
Dimension	150×120×115 (mm)
Net Weight	1.5 (kg)



#Let'sgive more 

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



Basic Power Supply

Series YR03410



Basic Power Supply

YR series Basic power supply comprises 2 models, including: YR03410 and YR03411, which can meet the conventional horizontal electrophoresis, small and medium-sized vertical electrophoresis, membrane acetate electrophoresis experiment requirement, the lightweight and compact, simple and practical characteristics, make it popular in for basic teaching.

YR03410 Basic Power Supply



Type of Output: Constant voltage or constant current
 Output range: JY300: 5~300V, 1~300mA, maximum 90W
 Increment: 1V, 1mA
 Timing Range: 1 min~9hr, 59min
 Display: LED
 Output jacks: Two sets of output jacks
 Dimension (LxWxH): 280x237x118 (mm)
 Weight: 3.0(kg)

- Molding shell, touch keys, microprocessor intelligent control;
- Automatic memory;
- Standard timing operation;
- Automatic detection of no-load, over-load.

YR03411 Basic Power Supplies



Type of Output: Constant voltage or constant current
 Output range: JY600: 5~600V, 1~200mA, maximum 120W
 Increment: 1V, 1mA
 Timing Range: 1 min~9hr, 59min
 Display: LED
 Output jacks: Two sets of output jacks
 Dimension (LxWxH): 280x237x118 (mm)
 Weight: 3.0(kg)

- Molding shell, touch keys, microprocessor intelligent control;
- Automatic memory;
- Standard timing operation;
- Automatic detection of no-load, over-load.



#Let'sgive more ♥


Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 Optimum Business Center 450 Rue Baden Powell,
 • 34000 Montpellier, France.
 Tit: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
A microscope is an instrument that allows observing and measuring objects that are too small to be seen with the naked eye. This is achieved through a complex system of lenses and mirrors, which form...

Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
A refrigerator is a device that is one of the most common pieces of laboratory equipment. It is used to maintain, in a controlled environment (refrigerated space), various fluids and substances...

Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



MONTPELLIER - FRANCE
HEADQUARTER
Optimum Business Center,
450 Rue Baden Powell,
34000 Montpellier - France.

Fax: +33 (0) 4 67 15 88 49
Tel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

Official Websites for Latin America.

