



LABORATORY SECTOR

PATHOLOGICAL ANATOMY



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Pathological Anatomy
Flotation Workstation (Water Bath-Slide Dryer)
Model YR456



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Pathological Anatomy Flotation Workstation (Water Bath-Slide Dryer)

Model YR456



Product Description:

- ✓ Pure-green digital screen for real-time monitoring of heating temperature, direct and clear, easy to operate; all parameters including preset temperatures, working temperatures, and wording status are real-time displayed
- ✓ Temperature is automatically program-controlled by single-chip microprocessors
- ✓ This device features multiple functions and its easy setting operation can meet the needs of different users;
- ✓ Two independent heating tank with DC low-voltage illuminating system and removable transparent heating dish: easy operation and convenient observation
- ✓ New heating mechanism using a new-type high-thermal-conductivity heating element provides even and quick heating: PID-controlled, Triple temperature controls, long lifespan, safe, reliable and energy-saving
- ✓ Temperature is precisely measured by temperature-sensing integrated blocks made in USA, and all settings are automatically stored in the system
- ✓ Special temperature-control and monitoring system characterized by strong resistance to abrasion and corrosion to ensure temperature precise;
- ✓ Automatic memory and restoration functions: After startup, all preset temperature data are automatically stored in the system
- ✓ Temperatures for floating and drying slides are individually controlled.

Technical Specifications:

Model	YR456
Left Temperature range of water bath	RT – 70°C; automatically maintained
Right Temperature range of water bath	RT – 70°C; automatically maintained
Temperature range of slide dryer	RT – 100°C; automatically maintained
Automatic memory and restoration functions	After startup, all preset temperature data
	are automatically stored in the system.
Temperature Control Precision	±2°C
Dimensions of water bath dish	155×155×55mm(W×D×H)
Dimensions of drying station(32pcs slides)	405×130mm(W×D×H)
Working Voltage	AC 220V±10% 50Hz (standard model); AC110V±10% 60Hz
Power	600W
Dimensions	475×460×175 mm (W×D×H)
Net weight	10kg



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

