



# LABORATORY LINES

## MOISTURE ANALYZERS



**Kalstein**®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



## KALSTEIN® Halogen Moisture Analyzer

Model Series YR05552



### Features :

- ✓ Halogen lamp/HBM sensor
- ✓ Aluminum shell/Stainless steel chamber
- ✓ LCD (White back light with black font)
- ✓ Temperature/Time set
- ✓ The percentage of moisture content/dry residual
- ✓ Store historical set
- ✓ Option: Aluminum box/Printer/Grinder



### Features YR053480A :

- ✓ Halogen lamp/HBM sensor
- ✓ Aluminum shell/Stainless steel chamber
- ✓ High-definition LCD
- ✓ Highlight halogen lamp
- ✓ Temperature/Time set
- ✓ The percentage of moisture content/dry residual
- ✓ Store historical set
- ✓ Option: Aluminum box/Printer/Grinder



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
 2 Rue Jean Lantier • 75001 Paris •  
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
 KALSTEIN FRANCE, S. A. S.



## KALSTEIN® Halogen Moisture Analyzer

Model Series YR05552



<b>Model</b>	YR05552	YR05553	YR05554	YR05555	YR05555A
<b>Capacity</b>			110g		
<b>Readability</b>	0.01g	0.005g	0.002g		0.001g
<b>Operating Temp</b>			5°C-35°C		
<b>Display</b>			LCD(White back light with black font)		LCD HD Back Display
<b>Calibration</b>			External Calibration		
<b>Interface</b>			RS232/RJ45/USB...		
<b>Pan Size</b>			Ø90mm		
<b>Packing Size</b>			490x350x360mm		
<b>Gross Weight</b>			7.5kg		
<b>Heating Source</b>			Halogen Lamp-1		Halogen Lamp-2
<b>Temp.Sensor</b>			PT100		
<b>Temp.Range</b>			40°C-199°C		
<b>Temp.Step</b>			1°C		
<b>Moisture Range</b>			0.00%-100.00%		
<b>Moisture Readability</b>	0.2%	0.10%	0.04%		0.01%
<b>Dry Range</b>			100.00%-0.00%		
<b>Dry Readability</b>	0.2%	0.10%	0.04%		0.01%
<b>Temp Setting</b>			40°C-199°C By 1°C Step		
<b>Time Setting</b>			1-99min By 10s		
<b>Number Of Storage</b>			15		
<b>Heater</b>			220V ±15% 50HZ/110V±15% 60HZ/400W		
<b>Balance</b>			Input:220V 50HZ/110V 60HZ Output:9V		



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S







 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación  
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale

<https://www.humatem.org/>



**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>



**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed

**The 3 articles of the week!**

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW

**How does the climate affect your microscope?**

microscope is an instrument that allows observing and measuring objects that are too small to be seen with the naked eye. This is achieved through a complex system of lenses and mirrors, which form

**Kalstein Blog**

**Article of the week:**  
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW

**What is the importance of the laboratory refrigerator?**

refrigerator is a device that is one of the most common pieces of equipment used in laboratories. It is a controlled environment (refrigerated space), various fluids and substances

**Kalstein Blog**

**Article of the week:**  
What is the importance of the laboratory refrigerator?


Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



**Kalstein France - Jacinto Convit Foundation / Alliances**

23 views • 3 years ago



**Kalstein's Microscopes**

445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**  
2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89  
**Cel:** +33 (0) 6 80 76 07 10  
[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

