




LABORATORY SECTOR

HEATING MANTLE / STIRRER



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Heating Mantle / Stirrer
Heating Mantle
Model YR05039 // YR05048



#Let'sgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Heating Mantle / Stirrer**Heating Mantle****Model YR05039 // YR05048****Product Description:**

YR heating mantle are the most widely used heating mantle in laboratories.

- ✓ Once shaping shell is made of cold rolled steel sheet with surface spraying treatment
- ✓ chrome heating wire, litter warming up time, long service life
- ✓ Easy operation, wide application
- ✓ Precise PID Temperature control meets high accuracy experiment needs
- ✓ Digital type single position mantles are with outer sensor, multi position mantles are with inner sensor, outer sensor optional.

**#Let'sgivemore**

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Heating Mantle / Stirrer**Heating Mantle****Model YR05039 // YR05048****Technical Description:**

Type	Analog				
Model	YR05039	YR05040	YR05041	YR05042	YR05043
Capacity (ml)	100	250	500	1 000	2 000
Voltage (V)	200 – 240V / 50 – 60Hz or 100 – 120V / 50 – 60Hz				
Max Use Temperature(°C)	380				
Speed (W)	100	150	250	350	450
Speed (r.p.m)	0 – 1 600rpm				
Working Time	Continu				
Exterior Size(mm)	Ø200×165		Ø260×200		Ø300×230
Packing Size(mm)	230 x 210 x 195		290x290x220		330x330x250
Net Weight (kg)	2.5		4.5		6.5
Gross Weight(kg)	3.5		5.5		7
Type	Digital				
Model	YR05044	YR05045	YR05046	YR05047	YR05048
Capacity (ml)	100	250	500	1 000	2 000
Voltage (V)	200 – 240V / 50 – 60Hz or 100 _ 120V / 50 – 60 Hz				
Max Use Temperature(°C)	380				
Speed (W)	100	150	250	350	450
Speed (r.p.m)	0 – 1 600rpm				
Working Time	Continu				
Exterior Size(mm)	Ø220 x 165		Ø260 x 200		Ø330 x 230
Packing Size(mm)	230x210x195		290x290x220		345x290x350
Net Weight (kg)	2.5		4.5		6.5
Gross Weight(kg)	3.5		5.5		7

**#Letsgivemore**

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:

 **Fundación
Jacinto Convit**
<https://www.jacintoconvit.org/>


 **ONE
TREE
PLANTED**

<https://onetreepanted.org>

 **Humatem**
Solidarité et coopération
biomédicale
<https://www.humatem.org/>


**MANIAPURE
FOUNDATION USA**
<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

1
Read this ARTICLE NOW



How does the climate affect your microscope?

microscope is an instrument that allows observing and measuring small objects. This is a very sensitive device that requires a controlled environment to maintain its performance.

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

2
Read this ARTICLE NOW



What is the importance of the laboratory refrigerator?

Refrigerator is a piece of equipment that is one of the most important in a laboratory. It is used to maintain samples in a controlled environment (refrigerated space), without frost and odors.

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

