



LABORATORY SECTOR

ANALYZER



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Analyzer
Hematology Analyzer For Veterinary
Model YR05122-V



#Let'sgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Analyzer

Hematology Analyzer For Veterinary

Model YR05122-V



Sample Profile:

- ✓ Dog, Cat, more species coming soon

Test Mode:

- ✓ CBC mode, CBC+DIFF mode
- ✓ Venous whole blood, Prediluted

Sample Volume:

- ✓ CBC+DIFF mode: ≤20ul
- ✓ CBC mode: ≤10ul

Reagents:

- ✓ Diluent (Vet)
- ✓ LH Lyse (Vet)
- ✓ LD Lyse (Vet)
- ✓ Probe Cleanser

Size / Weight:

- ✓ L*W*H = 480*375*517 (mm)
- ✓ Weight 36kg

Technical Specification:

- ✓ Tri-Angle Laser Scatter
- ✓ Flow Cytometer Method
- ✓ 3D Scattergram Analysis
- ✓ Impedance Method for RBC and PLT Counting
- ✓ Cyanide free reagent for HGB test

Parameters:

- ✓ 25 Reportable Parameters:
- ✓ WBC, RBC, HGB, HCT, MCV, MCH, MCHC, RDW-SD
- ✓ RDW-CV, PLT, MPV, PCT, PDW, P-LCR, P-LCC, NEU%
- ✓ LYM%, MON%, EOS,%, BAS%, NEU#, LYM#, MON#
- ✓ EOS#, BAS#
- ✓ 1 3D SCATTERGRAM
- ✓ 3 HIS HISTOGRAMS (WB C/ BASO, RBC, PLT)

Throughput:

- ✓ 60 test/hour

Data Memory:

- ✓ Up to 100,000 results (including histogram, Scattergram, animal information).

Display:

- ✓ 14 inch touch screen, resolution 1366*768

Interface:

- ✓ 1 LAN port, 4 USB ports

Printout:

- ✓ Support various external USB printers, printout formats user definable.

Communication:

- ✓ Support LIS connection, HL7 protocol formats user definable

Power Requirement:

- ✓ A.c.100-240V, 50/60Hz

Working Environment:

- ✓ Temperature 10-30 °C
- ✓ Humidity 20% - 85%
- ✓ Air Pressure: 70-106kPa
- ✓ Working Latitude: ≤350m



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



Fundación
Jacinto Convit

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale

<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
microscopes in a laboratory are one of the most sensitive pieces of equipment that require maintenance. This is because they are highly sensitive to temperature and humidity changes. This is why it is essential to maintain them in a controlled environment (refrigerated space), without dust and vibrations.

Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
Refrigerator in a laboratory is one of the most sensitive pieces of equipment that require maintenance. This is because they are highly sensitive to temperature and humidity changes. This is why it is essential to maintain them in a controlled environment (refrigerated space), without dust and vibrations.

Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

