




LABORATORY LINES

HOTPLATE MAGNETIC
STIRRER



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Laboratory Intelligent 5 - point Hotplate Magnetic Stirrer






Model YR02945



Features:

- ✓ Large Superfine black ceramic top plate 440 X 100 X 5MM for heating and stirring 5x 500ml at the same time
- ✓ Built-in PID electronic controller !
- ✓ Designed for stirring and heating of liquid. Built-in temperature control maintains element temperature from ambient to 320°C ! Stirring is Independently controlled with a variable speed control up to 1800rpm .
- ✓ LCD Display for Temperature set, reduce, add stir speed .
- ✓ Heat only, stir only or heat and stir together function.
- ✓ Electronic speed control for constant speed even during change in load .
- ✓ High magnetic adhesion .
- ✓ Self - tuning : You can set temperature as you want. It can stop working automatically when setting temperature has achieved. Temperature Control Accuracy $\pm 1^{\circ}\text{C}$.

Standard Accessories:

Photos	Item	Quantity
	Can be Removed Holder 1. Length 280mm 2. Diameter 10mm 3. Material : stainless steel .	1 PC
	Clamp with screw	1 PC
	Inner Thermocouple4 - spin	1 PC
	Outer Thermocouple4- spins Stainless Steel Length : 235mm Wire length : 1000mm	1 PC
	•PTFE-coated Length : 6mm Diameter : 20 mm	5 PCS



#Letsgivemore 

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.





All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Laboratory Intelligent 5 - point Hotplate Magnetic Stirrer

Model YR02945



Parameters

Model	YR02945
Temp Range for machine	RT ~V380°C
Temp Range for surface plate	RT~320 °C
Heat Dimension (L ×W × H) (mm)	440 X 100 X 5mm
Top Plate (L ×W × H) (mm)	470 X 140 X 5mm
Machine size	470x285x140mm
Electric Power (W)	600
Temperature Accuracy	± 1°C
Temperature Control	PID ,Inner and Outer Temperature Sensor
Plate Material	Superfine Black Ceramics
Stir Capacity for each point (mm)	5 X 50ml - 500ml
Motor Power for each point (W)	9W
Voltage	220V ,50Hz
Stir Speed	100 ~ 1800rpm
Stir display	± 1 rpm
Net Weight (g)	9.5 ± 0.2kg



#Let'sgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.




All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Laboratory Intelligent 5 - point Hotplate Magnetic Stirrer

Model YR02945



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:

 **Fundación
Jacinto Convit**
<https://www.jacintoconvit.org/>


 **ONE
TREE
PLANTED**

<https://onetreepanted.org>

 **Humatem**
Solidarité et coopération
biomédicale
<https://www.humatem.org/>

 **MANIAPURE
FOUNDATION USA**
<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
A microscope is an instrument that allows observing and measuring objects that are too small to be seen with the naked eye. This is achieved through a complex system of lenses and mirrors, which form...

Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
A refrigerator is a device that is one of the most common pieces of equipment used in laboratories. It is a controlled environment (refrigerated space), various fluids and substances...

Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

