


# Kalstein A different accompaniment, at your service

### **LABORATORY SECTOR**

MICROSCOPE





All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S





**OUR SERVICES** 

#### **Benefits and Support**

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



#### **Quick Response**

Our work team is always available to response all your consults or questions, in order to support you in any situation.





#### #Letsgivemore \*\*

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



#### **Technical Support**

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.





#### **Delivery Logistics**

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or.air.



#### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



#### **KALSTEIN® Microscope**

#### **Inverted Metallurgical Microscope**

**Model YR0258** 







#### #Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 /+33 6 80 76 07 10 • https://kalstein.eu

KALSTEIN FRANCE, S. A. S



#### **KALSTEIN® Microscope**

#### **Inverted Metallurgical Microscope**

#### **Model YR0258**






#### **Product Description:**

YR0258 inverted metallurgical microscope with stable support structure, advanced stage design and operation. It is widely used in institute and laboratory to identify the structure of various metals and other materials.

- ✓ Infinite optical system for high optical image quality
- ✓ Eyepiece Eye piece EW10x/20mm
- ✓ Infinity Plan Metallurgical Coverless achromatic objective
- ✓ Kohler 6V30W external lighting halogen lamp









#### #Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S. 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu





#### **KALSTEIN® Microscope**

#### **Inverted Metallurgical Microscope**

#### **Model YR0258**





#### **Technical Description**:

Model		YR0258		
Optical System	Infinity Optical System			
Head	Seidentopf Trinocular Head, 30° Inclined, Interpupillary 48~75mm, Diopter Adjustable			
Eyepiece	Extra Wide Field Eyepiece EW10x/20mm, High Eyepoint			
Nosepiece		Quintuple Nosepiece		
Objective	Infinity Plan MetallurgicalCoverless Achromatic Objective	4x/0.1/∞/-, W.D.17.3mm		
		10x/0.25/∞/-,W.D.10.2mm		
		20x/0.4/∞/-,W.D.5.1mm		
		40x/0.65/∞/-, W.D.0.52mm		
Focusing	Coaxial Coarse & Fine Adjustment, Coarse Stroke 37.7mm Per Rotation, Fine Stroke 0.2mm Per Rotation			
Working Stage	Plain Stage 160*250mm With Slide Clips			
Polarizing	Polarizing Attachment			
Filter	Green, Blue, Yellow, Ground Glass, Dia.45mm			
Illumination	External Kohler Illumination, Halogen Lamp 6V30W, Brightness Adjustable			
	Optional A	Accesories		
Objective	Infinity Plan MetallurgicalCoverless Achromatic Objective	5x/0.12/∞/-, W.D.15.4mm	YR0258-op1	
		50x/0.75/∞/-, W.D.0.32mm	YR0258-op2	
		80x/0.9/∞/-, W.D.0.20mm	YR0258-op3	
		100x/0.9/∞/0, W.D.0.7mm	YR0258-op4	
Wokring Stage		Attachable Mechanical Stage, X-Y Coaxial Control, Moving Range 120*78mm		
	Auxiliary Stage 70*180mm		YR0258-op6	
Adapter	Photo Adapter		YR0258-op7	
	C-Mount 1.0x		YR0258-op8	
	C-Mount 0.5x		YR0258-op9	



#Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 /+33 6 80 76 07 10 • https://kalstein.eu

KALSTEIN FRANCE, S. A. S







#LETSGIVEMORE >

# WITH THE ADCQUISITION OF A KALSTEIN EQUIPMENT

YOU MAKE A CONTRIBUTION TO:





https://onetreeplanted.org



Solidarité et coopération biomédicale

https://www.humatem.org/



https://www.maniapurefoundation.org/




## **ANY DOUBT?**Contact us!

#### **Articles Kalstein**

(Social Networks and Official Websites)







Article of the week: How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites https://kalstein.eu/articles





Article of the week: What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites https://kalstein.eu/articles



PARIS - FRANCE HEADQUARTER

2 Rue Jean Lantier 75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89 **Cel:** +33 (0) 6 80 76 07 10 sales@kalstein.eu

https://kalstein.eu/ https://www.kalstein.fr/

#### **Kalstein Training Videos**

(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago

