



LABORATORY SECTOR

MICROSCOPE



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Microscope
Laboratory Biological Microscope, Binocular, Achromatic
Model YR05787-1A



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S




KALSTEIN® Microscope
Laboratory Biological Microscope, Binocular, Achromatic
Model YR05787-1A



Product Description:

- ✓ Entry Level Laboratory Microscope With Full Range Upgrade Extension
- ✓ Can Use 40x LWD Infinity Plan Objective W.D.=3.0mm To Avoid Contaminated By Cedar Oil
- ✓ Creative 100x Water Objective Instead Of 100x Oil Objective For Easy Use Without Oil
- ✓ High Rigidity Body Structure Design Supply High Steady & Easy Handle
- ✓ Different Color Available To Customize Your Model



Seidentopf Head:

Inclined 30°, Rotatable 360°, Interpupillary Distance 48-75mm, High Eye point Eyepiece, Which Is Convenient For People With Different Sitting Height; The Integrated Design Of The Eyepiece Lens Barrel, Avoids The Loss And Contamination Of The Eyepiece And Enhances The Integrated Performance Of The Microscope.

Quadruple Backward:

4 Sets Objective Available: Achromatic, Plan, Infinity E-Plan, Infinity Plan.

Double Slide Holder Working Stage:

Mechanical Stage Size: 145 mm×140mm, Travel: 76mm×52mm, Scale: 0.1mm

Condenser and Light Source:

Abbe Condenser N.A.1.25, Iris Diaphragm, Rack & Pinion Up/Down, With Filter Holder 3W-LED Illumination Systems



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Microscope Laboratory Biological Microscope, Binocular, Achromatic Model YR05787-1A



Technical Specifications:

Model	YR05787-1A
Power Supply Type	DC USB Adapter
Transmit Light	3W LED, Brightness Adjustable
Light Source Type	LED
Condenser N.A.	1.25
Condenser Type	Abbe
Condenser Features	Rack & Pinion↕ Filter Holder Iris Diaphragm
Focusing Knob	Low Position Focusing Knob For Easy Operation
Focusing Scale	0.002mm
Focusing Range	26mm
Focusing Type	Coaxial Coarse & Fine
Focusing Features	Safety Stop Screw Tension Adjustable
Working Stage Moving	76x52mm
Working Stage Size	145x140mm
Working Stage Type	2 Layers Mechanical
Working Stage Features	Double Slide Holder Size ≥140x140mm
Objective	Achromatic Objective 4x,10x,40x,100x
Objective Type	Achromatic
Nosepiece Holes	4 Holes
Nosepiece Features	Inward
Eyepiece	WF10x/18mm, 2 pcs
Eyepiece Diameter	23.2mm
Head Diopter	Left Adjustable
Head Interpupillary	48~75mm
Head Rotation	Head Rotate 360°
Head Inclination	30°
Standard Magnification	40x-1000x
Optical System	Achromatic
Upgrade Function	BF+PH+DF+PL+FL
Standard Function	BF



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

