


# Kalstein A different accompaniment, at your service

## **LABORATORY SECTOR**

**BIOSECURITY CABINETS** 





All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 /+33 6 80 76 07 10 • https://kalstein.eu

KALSTEIN FRANCE, S. A. S





**OUR SERVICES** 

### **Benefits and Support**

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



# Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



#### **Quick Response**

Our work team is always available to response all your consults or questions, in order to support you in any situation.





#### #Letsgivemore \*\*

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



#### **Technical Support**

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.





#### **Delivery Logistics**

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or.air.



#### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



### **KALSTEIN® Biosecurity Cabinets**

**Laminar Flow Cabinet** 

**Model YR05738** 







#### #Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S. 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu

KALSTEIN FRANCE, S. A. S



#### **KALSTEIN® Biosecurity Cabinets**

**Laminar Flow Cabinet** 

**Model YR05738** 





#### **Product Description:**

This is a kind of clean bench widely used in medical treatment, chemical experiment, electron, precise machinery, instrument, apparatus, medical preparation, operating room, foodstuff, drinks, daily use cosmetic, stock raising and aviation, as well as providing germfree dustless clean environment.

- ✓ Adopting arbitrary positioning sliding door system
- ✓ The shell is welded by cold-rolled plates, the surface is electrostatically sprayed, and the work surface is SUS304 brushed stainless steel, which is corrosion-resistant and easy to clean.
- ✓ The lighting and sterilization systems are interlocked safely, and have the timing function of the ultraviolet lamp
- ✓ Digital display LCD contrail interface, more user-friendly design.
- ✓ The vertical quasi-closed table top and the formation of a downward flow air curtain in the operation room con effectively prevent the input of external air and the operation area is clean.
- Equipped with a HEPA filter and a primary filter for preliminary filtration, which con effectively extend the service lite of the high-efficiency filter.

#### **Technical Specifications:**

Model	YR05738
Air supply	Vertical Air Suppply
Cleanliness grade	Class 100 (US Federal 209E)
Number of colonies	≤0.5pcs/dish-hour (Ø90mm culture dish)
The average wind speed	0.25-0.45m/s
Noise	≤62dB (A)
Power supply	AC single phase 220V /50Hz
Vibration half peak	≤0.5µM (x, y, z direction)
Illumination	≥300LX
Max. Power consumption	400 W
Dimension of working area	000.050.500
(W*D*H)(mm)	880x650x580
Overall dimension (W+D+H)(mm)	1040x690x1600
Specification and quantity	
of high efficiency filter	865x555x50x 1
Specification and quantity	
	12 W x 1 / 30 W x 1
of fluorescent lamp/UV lamp	
Applicable persons	Single person- Single Side



#### #Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu

KALSTEIN FRANCE, S. A. S







#LETSGIVEMORE >

# WITH THE ADCQUISITION OF A KALSTEIN EQUIPMENT

YOU MAKE A CONTRIBUTION TO:





https://onetreeplanted.org



Solidarité et coopération biomédicale

https://www.humatem.org/



https://www.maniapurefoundation.org/



# **ANY DOUBT?**Contact us!

#### **Articles Kalstein**

(Social Networks and Official Websites)







Article of the week: How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites https://kalstein.eu/articles





Article of the week: What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites https://kalstein.eu/articles



PARIS - FRANCE HEADQUARTER

2 Rue Jean Lantier 75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89 **Cel:** +33 (0) 6 80 76 07 10 sales@kalstein.eu

https://kalstein.eu/ https://www.kalstein.fr/

#### **Kalstein Training Videos**

(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago

