




# LABORATORY SECTOR

PIPETTES



## Kalstein®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



**KALSTEIN® Pipettes**  
**Liquid Handling Pipettes**  
**Model YR03733 // YR03733-8**



**#Letsgivingore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
**KALSTEIN FRANCE, S. A. S**



## KALSTEIN® Pipettes Liquid Handling Pipettes Model YR03733 // YR03733-8




### Product Description:

YR Pipettes Series are a new mechanical, air-displacement pipette that offers you a volume lock and this one of the most beautiful exterior in pipettes. Enjoy a new addition to pipette line with its legendary accuracy, precision, and robustness... with the highest comfort.

### Features:

- ✓ 3 digits
- ✓ Magnifying Shape
- ✓ Chemical resistance
- ✓ 121°C Sterilization
- ✓ Very low operating force
- ✓ Color indicates pipette volume
- ✓ Positioned for perfect ergonomics
- ✓ Only a few turns to reach your desired volume



Model	Volume Range	Accuracy	Repeatability
YR03733	0.1 – 2 µl	±12.0 – 1.5%	≤6.00 – 0.70%
YR03733-1	0.5 – 10 µl	±2.5 – 1.0 %	≤1.50 – 0.40%
YR03733-2	2 – 20 µl	±2.5 – 1.0 %	≤1.50 – 0.30%
YR03733-3	5 – 50 µl	±2.0 – 1.0%	≤1.50 – 0.25%
YR03733-4	10 – 100 µl	± 1.8 – 0.8%	≤0.50 – 0.15%
YR03733-5	20 – 200 µl	±1.8 – 0.6%	≤0.50 – 0.15%
YR03733-6	100 – 1 000 µl	±1.5 – 0.6%	≤0.30 – 0.15%
YR03733-7	500 – 5 000 µl	±1.2 – 0.6%	≤0.30 – 0.15%
YR03733-8	1 – 10ml	±2.5 – 0.6%	≤0.30 – 0.15%



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S







 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación  
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale

<https://www.humatem.org/>



**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>



**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed  
**The 3 articles of the week!**  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW  
**How does the climate affect your microscope?**  
Microscopes are instruments used for observing microscopic objects. This is a sensitive device that requires a controlled environment to maintain its performance.  
**Kalstein Blog**  
**Article of the week:**  
How does the climate affect your microscope?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW  
**What is the importance of the laboratory refrigerator?**  
Refrigeration is a key element of the most advanced equipment that is used in laboratories. It is a controlled environment (refrigerated space), without fluctuations in temperature.  
**Kalstein Blog**  
**Article of the week:**  
What is the importance of the laboratory refrigerator?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



**Kalstein France - Jacinto Convit Foundation / Alliances**  
23 views • 3 years ago



**Kalstein's Microscopes**  
445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**  
2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89  
**Cel:** +33 (0) 6 80 76 07 10  
[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

