




MEDICAL LINES MESH NEBULIZER



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.





Mesh Nebulizer

Model YR21A



Characteristics:

- ✓ Microporous direct atomization technology, tiny particles(3-5 microns)
- ✓ Ultra-quiet, low power consumption, built-in lithium battery
- ✓ Almost without residual liquid
- ✓ Mini design, easy to carry
- ✓ Two modes of automatic turn off by 5 or 10 minutes

Package Included:

- 1 x Ultrasonic Mesh Nebulizer
- 1 x Adult Mask
- 1 x Child Mask
- 2 x Mouthpiece
- 1 x USB charging line
- or
- 1 x Adult Mask
- 1 x Child Mask
- 2 x Mouthpiece

Technical Specifications:

- ✓ Product Name: Ultrasonic Mesh Nebulizer
- ✓ Material: ABS
- ✓ Size: 42x55x119mm
- ✓ Model: YR21A
- ✓ Atomization Rate: $\geq 0.3\text{ml/min}$
- ✓ Power: DC 3.7V lithium battery
- ✓ Battery Capacity: 800mAh
- ✓ Power: Approx.1.5W
- ✓ Particle Size: $2.7\mu\text{m}\pm 25\%$
- ✓ Medicine Temperature: $\leq 60^\circ\text{C}$
- ✓ Ultrasonic Vibration Frequency: 110kHz
- ✓ Noise Of the Whole Machine: $\leq 50\text{dB (a)}$
- ✓ Medicine Capacity: 20ml maximum,0.2ml minimum
- ✓ Working Environment Temperature: $5^\circ\text{C} \sim 40^\circ\text{C}$
- ✓ Relative Humidity: $15\% \sim 90\% \text{RH}$, not in Coagulation situation Air Pressure: (70.0 ~ 106.0)kPa
- ✓ Battery Life: the battery can be used for approximately 5 days after charging (if operating the device for 30 minutes a day and the battery cycle life is 300 cycles)



#Letsgivemore

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S





Mesh Nebulizer

Model YR21A

Features:

- ✓ The nebulizer is a device specially used to treat upper respiratory diseases, such as cold,
- ✓ fever, cough, phlegm, asthma, sore throat, pharyngitis, rhinitis, bronchitis, pneumonia
- ✓ and other diseases.
- ✓ The drug enters the respiratory tract and the lungs by breathing inhalation and reaches the affected part, so the atomization treatment works faster than the oral medication.
- ✓ The atomization treatment has a quick effect and high safety.
- ✓ Microporous direct atomization technology
- ✓ Ultra-quiet, low power consumption
- ✓ Built-in lithium battery
- ✓ Almost without residual liquid
- ✓ Two modes of automatic turn off by 5 or 10 minutes



Nebulizer

Bite

USB charging cable

Adult face mask

Child mask



#Letsgivemore ♥


Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 Optimum Business Center 450 Rue Baden Powell,
 • 34000 Montpellier, France.
 Titl: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
A microscope is an instrument that allows observing and measuring objects that are too small to be seen with the naked eye. This is achieved through a complex system of lenses and mirrors, which form...

Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
A refrigerator in a laboratory is one of the most important equipment. It is used to maintain, in a controlled environment (refrigerated space), various fluids and substances...


Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



MONTEPELLIER - FRANCE
SEDE

Optimum Business Center,
450 Rue Baden Powell,
34000 Montpellier - France.

Fax: +33 (0) 4 67 15 88 49

Tel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

Official Websites for Latin America.

