


# Kalstein A different accompaniment, at your service

# LABORATORY LINES MICROPIPETTES





All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu

KALSTEIN FRANCE, S. A. S





OUR SERVICES

### **Benefits and Support**

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### **Quick Response**

Our work team is always available to response all your consults or questions, in order to support you in any situation.





### #Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breacancer and helps most vulnerable communities.



### **Technical Support**

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.





### **Delivery Logistics**

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or.air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.





### **Single Channel Adjustable Volume Pipette**



### Series YR03735

Kalstein's pipette series YR03757, is the latest innovation of Kalstein based on ergonomics, comfort, and traceability research. This new mechanical pipette combines the Kalstein's pipette accuracy, reliability and newly designed body.

### Features:


- ✓ 4 Digits
- ✓ Magnifying shape
- ✓ Chemical Resistance
- ✓ 121°C sterilization
- ✓ Very low operation force
- ✓ Color indicates pipette volume
- ✓ Positioned for perfect ergonomics
- ✓ Only a few turns to reach your desired volume

#### Remark:

- ✓ Quantity per Carton: 20pcs
- ✓ G.W: 7,0 KG
- ✓ Dimension: 43 x 31 x 27 (cm)

### Control button

- · Very low operating force
- · Colour indicates pipette volume





### #Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S. 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu

KALSTEIN FRANCE, S. A. S





# Single Channel Adjustable Volume Pipette



### Series YR03735

Kalstein's pipette series YR03757, is the latest innovation of Kalstein based on ergonomics, comfort, and traceability research. This new mechanical pipette combines the Kalstein's pipette accuracy, reliability and newly designed body.

#### Features:

- ✓ 4 Digits
- ✓ Magnifying shape
- ✓ Chemical Resistance
- ✓ 121°C sterilization
- ✓ Very low operation force
- ✓ Color indicates pipette volume
- ✓ Positioned for perfect ergonomics
- ✓ Only a few turns to reach your desired volume

### Remark:

✓ Quantity per Carton: 20pcs

✓ G.W: 7,0 KG

 $\checkmark$  Dimension: 43 x 31 x 27 (cm)

Model	Volume Range	Accuracy	Repeatability
YR03735	0.1-2µ	± 12.0-1.5%	≤6.00-0.70%
YR03735-1	0.5-1 OµI	± 2.5-1.0%	≤1.50-0.40%
YR03735-2	2-20μ1	± 2.5-1.0%	≤1.50-0.30%
YR03735-3	5-50μ1	± 2.0-1.0%	≤1.50-0.25%
YR03735-4	10-100µ1	± 1.8-0.8%	≤0.50-0.15%
YR03735-5	20-200µ1	± 1.8-0.6%	≤0.50-0.15%
YR03735-6	100-1000μΙ	± 1.5-0.6%	≤0.30-0.15%
YR03735-7	500-5000µ1	± 1.2-0.6%	≤0.30-0.15%
YR03735-8	1-10ml	± 2.5-0.6%	≤0.30-0.15%



#Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu

KALSTEIN FRANCE, S. A. S





## Multichannel Adjustable Volume Pipette



### Series YR03736

Kalstein's pipette series YR03764 and YR03765, is the latest innovation of Kalstein based on ergonomics, comfort, and traceability research. This new mechanical pipette combines the Kalstein's pipette accuracy, reliability and newly designed body.

### **Product** Features

- ✓ 3 Digits
- ✓ Magnifying shape
- ✓ Chemical Resistance
- √ 121°C sterilization
- ✓ Very low operation force
- ✓ Color indicates pipette volume
- ✓ Positioned for perfect ergonomics
- ✓ Only a few turns to reach your desired volume

### Remark

✓ Quantity per Carton: 10pcs

✓ G.W: 8,0 KG

✓ Dimension: 69 x 32 x 23 (cm)

Model	Channels	Volume Range	Accuracy	Repeatability
YR03736	8	2-20µ	± 6.0-1.5%	<b>≤</b> 2.50-0.40%
YR03736-1	8	5-50µ	± 3.0-1.2%	≤1.50-0.40%
YR03736-2	8	20-200μ	± 2.5-1.0%	≤1.25-0.50%
YR03737	12	2-20µ	± 6.0-1.5%	<b>≤</b> 2.50-0.40%
YR03737-1	12	5-50µ	± 3.0-1.2%	≤1.50-0.40%
YR03737-2	12	20-200µ	± 2.5-1.0%	≤1.25-0.50%



### #Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu

KALSTEIN FRANCE, S. A. S







#LETSGIVEMORE >

# WITH THE ADCQUISITION OF A KALSTEIN EQUIPMENT

YOU MAKE A CONTRIBUTION TO:





https://onetreeplanted.org



biomédicale

https://www.humatem.org/



https://www.maniapurefoundation.org/



## **ANY DOUBT?**Contact us!

### **Articles Kalstein**

(Social Networks and Official Websites)







Article of the week: How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites





Article of the week: What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites https://kalstein.eu/articles



PARIS - FRANCE HEADQUARTER

2 Rue Jean Lantier 75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89 **Cel:** +33 (0) 6 80 76 07 10 sales@kalstein.eu

https://kalstein.eu/ https://www.kalstein.fr/


### **Kalstein Training Videos**

(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago

