




# LABORATORY LINES

## HOTPLATES



**Kalstein**®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



# KALSTEIN® HOTPLATES

## Multi-position magnetic stirrer

Model YR02955



### Features

1. 4-position Hotplate Magnetic stirrer
2. Digital Display, Temperature: Rt -380°C, Mixing speed :100-1800 RPM
3. Ceramic-coated aluminum plate for corrosion resistance
4. Specially designed housing prevents spill from damaging the electronic components inside
5. Adjustable rotation speed and temperature

### Specification

Model	YR02955
Top Plate Dimension (mm)	140x 140mm (5.5x5.5inch )
Temperature Range for Inner Plate (°C )	400 degree
Temperature Range for Plate Surface (°C )	330 degree
Position	6
Heat Control	PID
Thermocouple type	PT100
Thermocouple length	SS316 :170mm.Wire:1000mm
Top Material	Black ceramics, shape : round or square plate
Temperature Display (°C )	+/-1
Stir Speed	100 ~ 1500rpm
Stir Display	Digital Display
Stir speed Accuracy	+/- 1 rpm
Permissible Ambient Temperature	Ambient Temperature to 40°C
Shell Material	Superfine Aluminum alloy
Voltage	220v ,50Hz / 110v 60 Hz
Heating Power(w)	2700W
Motor Power for Each Position	40W
Machine Power	2940W
Temperature Accuracy (°C )	+/-1
Stirring Capacity of Each Position (Water)	50- 2000ml



#Letsgivingmore 

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
Optimum Business Center 450 Rue Baden Powell,  
• 34000 Montpellier, France.  
Tit: +33 467158849 / +33 680760710/ +33 663810023  
<https://kalstein.eu>  
KALSTEIN FRANCE, S.A. S







 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:

 **Fundación  
Jacinto Convit**  
<https://www.jacintoconvit.org/>

 **ONE  
TREE  
PLANTED**


<https://onetreepanted.org>

 **Humatem**  
Solidarité et coopération  
biomédicale  
<https://www.humatem.org/>

  
**MANIAPURE  
FOUNDATION USA**  
<https://www.maniapurefoundation.org/>



**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed  
**The 3 articles of the week!**  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW  
**How does the climate affect your microscope?**  
microscopes in a laboratory are one of the most sensitive pieces of equipment that require to maintain in a controlled environment (refrigerated space), without fluctuations in humidity and temperature.  
**Kalstein Blog**  
**Article of the week:**  
How does the climate affect your microscope?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW  
**What is the importance of the laboratory refrigerator?**  
Refrigerator in a laboratory is one of the most sensitive pieces of equipment that require to maintain in a controlled environment (refrigerated space), without fluctuations in humidity and temperature.  
**Kalstein Blog**  
**Article of the week:**  
What is the importance of the laboratory refrigerator?  
Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



**Kalstein France - Jacinto Convit Foundation / Alliances**  
23 views • 3 years ago



**Kalstein's Microscopes**  
445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**  
2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89  
**Cel:** +33 (0) 6 80 76 07 10  
[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

