



MEDICAL SECTOR

OPERATION TABLES



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Operation Tables
Operation Tables
Model YR02154



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Operation Tables**Operation Tables****Model YR02154****Features:**

- ✓ YR02154 surgical operating table features X-ray radiolucent surface and allows convenient use of C-arm.
- ✓ It uses conformable pads that are easy to replace and clean.
- ✓ The surgical operating table is highly maneuverable with a central brake.
- ✓ It adopts the integrated rechargeable battery that allows 50 nonstop batches of surgical operations after one full charge.
- ✓ Longitudinal shift of up to 296 mm is available

Product Description:

YR02154 surgical operating table is a universal electric operating table that allows all kinds of surgical positions. Based on the highly reliable and stable electro-hydraulic system, it can be easily set up and precisely positioned with ultra low noise. With integrated one key reset function, handheld panel, standby control panel in the column and foldable head plate for greater rest, the surgical operating table will bring great operation convenience and surgical experience for both surgeons and patients.

Application:

YR02154 is an electro-hydraulic operating table designed to meet universal surgical-positioning needs, like in Neurology, Ophthalmology, Gynecology, Urology, Minimal Invasive Surgery, Hand & Arm surgery .

Model	YR02154
Table length	2100 mm
Table width	500mm (not include side rail)
Table height (up / down)	1030mm / 700mm
Trendelenburg / Reverse Trendelenburg	30 ° / 30 °
Lateral tilt (right / left)	20 ° / 20 °
Back plate (down)	68 ° / 42 °
Head plate (down)	57 ° / 72 °
Leg plate (downsplit)	96 ° / 65 ° x 2
Flexion / Reflection	220 ° / 115 °
Longitudinal displacement	296 millimeters
Patient weight	200Kg (during positioning)
Input power	AC 110 / 220V, 60 / 50Hz

**#Letsgivemore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

