



LABORATORY SECTOR

PH METERS



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® PH Meters
Pocket Salinity Meter
Model YR01815-1



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® PH Meters
Pocket Salinity Meter
Model YR01815-1



Product Description:

High-range pocket salinity tester, 1 to 3 points calibration, setup menu contains 5 options. The meter is suitable for measuring the seawater, accuracy: 1% F.S.

- ✓ High-range pocket salinity tester contains the conductivity measurement mode.
- ✓ Platinum conductivity cell provides the quick and reliable measurement results.
- ✓ 1 to 3 points calibration with automatic recognition for conductivity standards.
- ✓ Automatic temperature compensation ensures accurate readings over the entire range.
- ✓ Auto-Read function senses and locks the measurement endpoint.
- ✓ Auto-Power Off effectively conserves battery life.
- ✓ Setup menu allows setting the number of calibration points, temperature unit, etc.
- ✓ Reset function automatically resumes all settings back to the factory defaults.
- ✓ Replaceable electrode module reduces maintenance and replacement cost.
- ✓ Waterproof feature to ensure complete protection in harsh environments.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® PH Meters
Pocket Salinity Meter
Model YR01815



Product Description:

Low-range pocket salinity tester, 1 to 3 points calibration, setup menu contains 5 options. The meter is suitable for general applications, accuracy: 1% F.S.

- ✓ Low-range pocket salinity tester contains the conductivity measurement mode.
- ✓ Platinum conductivity cell provides the quick and reliable measurement results.
- ✓ 1 to 3 points calibration with automatic recognition for conductivity standards.
- ✓ Automatic temperature compensation ensures accurate readings over the entire range.
- ✓ Auto-Read function senses and locks the measurement endpoint.
- ✓ Auto-Power Off effectively conserves battery life.
- ✓ Setup menu allows setting the number of calibration points, temperature unit, etc.
- ✓ Reset function automatically resumes all settings back to the factory defaults.
- ✓ Replaceable electrode module reduces maintenance and replacement cost.
- ✓ Waterproof feature to ensure complete protection in harsh environments.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® PH Meters

Pocket Salinity Meter

Model YR01815



Technical Description

Model	YR01815
Salinity	
Range	0.00 to 10.00ppt
Resolution	0.01ppt
Accuracy	±1% F.S.
Conductivity	
Range	0 to 20.00, 200.0, 2000µS/cm, 20.00mS/cm
Resolution	0.01, 0.1, 1
Accuracy	±1% F.S.
Calibration Points	1 to 3 points
Calibration Solutions	84µS, 1413µS, 12.88mS
Temperature	
Range	0 to 60°C/32 to 140°F
Resolution	0.1°C/0.1°F
Accuracy	±1°C/±1.8°F
Offset Calibration	1 point
Calibration Range	Reading ±10°C
General Specifications	
Temperature Compensation	0 to 60°C/32 to 140°F, automatic
Temperature Coefficient	2%/°C, linear
Reference Temperature	25°C
Cell Constant	K=1
Hold Function	Manual or auto-endpoint
Auto-Power Off	8 minutes after last key pressed
Sensor Type	ECSCAN-C1-10K
Display	Dual-line LCD (21×21mm)
Power	2×1.5V AAA batteries
Battery Life	Approximately 200 hours of continuous use
Dimensions	185(L)×40(Dia.)mm
Weight	100g



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

