




# LABORATORY SECTOR

MICROSCOPE



## Kalstein®

A different accompaniment, at your service





# Kalstein®

A different accompaniment, at your service

OUR SERVICES

## Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



### Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



### Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



### #Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



### Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



### Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



### Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



## KALSTEIN® Microscope Polarizing Microscope Model YR0264 // YR0265



**#Letsgivemore** ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S



## KALSTEIN® Microscope Polarizing Microscope Model YR0264 // YR0265



### Product Description:

- ✓ Polarizing microscope for normal inspection and education.
- ✓ Widely used in crystallography, Geology, Analytical chemistry, Medicine, Biology, Environmental ology, Pharmacy, Toxicology, Asbestos analysis and so on.
- ✓ Single polarization, Orthogonal polarization or horoscopic observation is available.



#Letsgivingore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S



## KALSTEIN® Microscope Polarizing Microscope Model YR0264 // YR0265



### Technical Specifications:

Model		YR0264	YR0265
<b>Head</b>	Sliding Binocular Viewing Head, Inclined at 45°	•	–
	Sliding Trinocular Viewing Head, Inclined at 45°	–	•
<b>Eyepiece</b>	Eyepiece WF10×/18 with scale of cross hair	•	•
	Eyepiece WF10×/18 with graticule	•	•
	WF10×/18	•	•
<b>Objective</b>	Strain Free Achromatic Objective 4x, 10x, 25x, 40x, 63x	•	•
	Strain Free Plan Achromatic Objective 25x, 40x, 63x	•	•
<b>Nosepiece</b>	Backward quadruple nosepiece, Center adjustable	•	•
<b>Focusing</b>	Coaxial coarse & fine adjustment, Range 28mm, Fine division 0.002mm	•	•
<b>Stage</b>	Revolving Round Stage Diameter $\Phi$ 160mm, Graduated in 1° Increments, Minimum resolution 6' when using vernier scale	•	•
	Mechanical stage, Moving range 30X40mm	•	•
<b>Condenser</b>	NA1.25 Abbe condenser with iris diaphragm & filter	•	•
<b>Analyzer</b>	Rotatable analyzer with gradation 0°-90°	•	•
<b>Polarizer</b>	360° Rotatable	•	•
<b>Bertrand Lens</b>	Bertrand Lens, sliding in/out of optical path	•	•
<b>Optical Compensator</b>	$\lambda$ Slip (First class red) , 1/ 4 $\lambda$ Slip, Quartz wedge	•	•
<b>Illumination</b>	12V/20W Halogen lamp, Brightness adjustable (Transmitted illumination)	•	•
	12V/50W Halogen lamp, Brightness adjustable (Reflected illumination)	•	•
<b>Micrometer</b>		•	•



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,  
2 Rue Jean Lantier • 75001 Paris •  
+33 1 78 95 87 89 / +33 6 80 76 07 10 •  
<https://kalstein.eu>  
KALSTEIN FRANCE, S. A. S





 **Kalstein**<sup>®</sup>  
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE  
ACQUISITION  
OF A KALSTEIN  
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación  
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



**Humatem**

Solidarité et coopération  
biomédicale


<https://www.humatem.org/>



**MANIAPURE  
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

**Articles Kalstein**  
(Social Networks and Official Websites)



**Kalstein Blog**  
Not just be informed, be fully informed

**The 3 articles of the week!**

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**1**  
Read this ARTICLE NOW

**How does the climate affect your microscope?**

**Kalstein Blog**

**Article of the week:**  
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>



**2**  
Read this ARTICLE NOW

**What is the importance of the laboratory refrigerator?**

**Kalstein Blog**

**Article of the week:**  
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites  
<https://kalstein.eu/articles>

**Kalstein Training Videos**  
(Youtube Channel and Official Websites)



**Kalstein France - Jacinto Convit Foundation / Allianc...**  
23 views • 3 years ago



**Kalstein's Microscopes**  
445 views • 1 week ago



PARIS - FRANCE  
**HEADQUARTER**  
2 Rue Jean Lantier  
75001 Paris - France

**Tlf:** +33 (0) 01 78 95 87 89  
**Cel:** +33 (0) 6 80 76 07 10  
[sales@kalstein.eu](mailto:sales@kalstein.eu)

<https://kalstein.eu/>  
<https://www.kalstein.fr/>

