




LABORATORY LINES

LAMINAR FLOW CABINET



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Laminar Flow Cabinet Single Person Laminar Flow Cabinet Model Series YR05266



Product Description

Can choose horizontal and vertical flow, single person single side; Mirror stainless steel, corrosion-resistant; low consumption, high efficient filter; microcomputer control, weight-balance front door;

Is the essential equipment to carry out biotechnology research and experimental, widely used in medical, college and other industries. provide air purification for the working area.

Standard Features

- ✓ Two models horizontal laminar flow and vertical laminar clean bench.
- ✓ According to using function are desktop and vertical clean bench, single person single side (two sides)
- ✓ The working surface is made of mirror stainless steel, easy to clean, corrosion-resistant.
- ✓ Equipped with many kinds filter to meet the customer's purification degree requirements.
- ✓ LED microcomputer controller, adjustable fan speed.
- ✓ Operation door can be located any position



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Laminar Flow Cabinet
Single Person Laminar Flow Cabinet
Model Series YR05266



Technical Parameters

Model	YR05266
Wind speed (m/s)	0.25 ~ 0.45
Illumination (LX) & Noise	Illumination ≥300lx
Voltage (V)	220
Power (KW)	0.3
Inner size (mm) (H×W×D)	500×700×430
Overroll size (mm) (H×W×D)	1610×860×580
Packing size (mm) (H×W×D)	1200×940×660
Net weight (KG)	80
Shipping weight (KG)	85



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Single Person Laminar Flow Cabinet Single Person Laminar Flow Cabinet Model Series YR05266



Catalog of Biosafety Campaigns and Booths on offer



Stainless Steel Laboratory Fume Hood YR05E

Select Options



Laminar Flow Cabinet YR05746

Add To Quote List



Epoxy-Coated Steel Basic Chemical Laboratory

Select Options



Laboratory Fume Hood, acid and alkali resistant

Select Options



Horizontal Laminar Flow Cabinet YR05735

Add To Quote List



Laminar Flow Cabinet YR05737

Add To Quote List



Biological Safety Cabinet 100% Exhaust Class

Select Options



Laminar Flow Cabinet YR05743

Add To Quote List



#Letsgivingore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



Fundación
Jacinto Convit

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

1
Read this ARTICLE NOW



How does the climate affect your microscope?

microscope is an instrument that allows observing and measuring small objects. This is possible through the use of lenses or a system of mirrors that collect and focus light from the object being observed.

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

2
Read this ARTICLE NOW



What is the importance of the laboratory refrigerator?

Refrigerator is a device that is one of the most common equipment that is used to maintain, in a controlled environment (refrigerated space), various fluids and substances.

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

