


Kalstein A different accompaniment, at your service

LABORATORY SECTOR

REFRIGERATORS AND FREEZERS





All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 /+33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S





OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.





#Letsgivemore **

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.





Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or.air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Refrigerators and Freezers

Top Open 305L Laboratory Deep Freezer

Model YR05317







#Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 /+33 6 80 76 07 10 • https://kalstein.eu

KALSTEIN FRANCE, S. A. S



KALSTEIN® Refrigerators and Freezers

Top Open 305L Laboratory Deep Freezer

Model YR05317





Product Description:

Laboratory Deep Freezer Top Open With Microcomputer Control For Lab And Hospital Use

Promed Deep Freezer are designed to store samples products.

- ✓ Wide temperature range :10°C~32°Cambient temperature can reach -10°C ~-25°C
- ✓ Accurate Temperature control technology to make the temperature maintain-10°C ~-25°Cand accuracy 0.1 °C.
- ✓ SPCC sprayed steel external material and inner chamber for friendly using and easy cleaning.
- ✓ High-density polyurethane foam(65mm)

Safety

- Build-in back-up battery to display temperature and to operate audible and visual alarm system for up to 12 hours without AC power.
- ✓ Visual and audible alarm system: High/low temperature, Sensor error, Thermostat failure alarm, power failure alarm
- ✓ Power failure protection: turn on delay of the cooling system when power failure, restart delay protection.
- Keyboard lock, password protection to avoid set parameter randomly.
- ✓ Wide available voltage range: 198V~242V.

Humanized Design:

- ✓ 4 casters,2 casters with stopper for easy moving and fixation.
- ✓ 1 Standard test hole(25mm),convenient for temperature monitoring.

Certificate:

- ✓ ISO14001.ISO9001.ISO13485.ISO45001
- ✓ FDA List and CE

Warranty:

✓ 1 year free spare parts and remote technical support.

Special:

✓ 4 visual and audible alarms. 1 Standard test hole(25mm).







#Letsgivemore 💗

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S. 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 / +33 6 80 76 07 10 • https://kalstein.eu KALSTEIN FRANCE, S. A. S



KALSTEIN® Refrigerators and Freezers

Top Open 305L Laboratory Deep Freezer









Technical Description:

Model	YR05317
	Performance
Temperature Range	-10~-25°C
Ambient Temperature	+10~+32°C
Controller	Microprocessor
Display	Digital Display
R	Refrigeration System
Compressor/Brand	1/China
Refrigerant	R134A
Cooling Type	Direct cooling
Defrost Type	Manual
	Electrical Data
Power Supply(V/Hz)	220~240V/50~60;115/60
Rated Power(W)	106
Electrical Current(A)	0.92
Power Consumption	1.5 kWh/24h
	Dimensions
Capacity	305L
Interior Size(W*D*H)	1000*520*665mm
Exterior Size(W*D*H)	1135*725*890mm
Package Size(W*D*H)	1223*777*1005mm
N.t./G.t.(approx)	74.2/87kg
Loading Qty(20'/40'/40'H)	26/56/56
	Alarm
High/Low Temperature	Υ
Sensor Error	Υ
Power Failure	Y
	Accessories
Caster	4
Remote alarm port	1
Test hole	1/25MM





Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved ® KALSTEIN France S. A. S., 2 Rue Jean Lantier • 75001 Paris • +33 1 78 95 87 89 /+33 6 80 76 07 10 • https://kalstein.eu



KALSTEIN FRANCE, S. A. S





#LETSGIVEMORE >

WITH THE ADCQUISITION OF A KALSTEIN EQUIPMENT

YOU MAKE A CONTRIBUTION TO:





https://onetreeplanted.org



Solidarité et coopération biomédicale

https://www.humatem.org/



https://www.maniapurefoundation.org/



ANY DOUBT?Contact us!

Articles Kalstein

(Social Networks and Official Websites)







Article of the week: How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites https://kalstein.eu/articles





Article of the week: What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites https://kalstein.eu/articles



PARIS - FRANCE HEADQUARTER

2 Rue Jean Lantier 75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89 **Cel:** +33 (0) 6 80 76 07 10 sales@kalstein.eu

https://kalstein.eu/ https://www.kalstein.fr/

Kalstein Training Videos

(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago

