




LABORATORY LINES

ULTRASONIC BATHS



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Ultrasonic Baths

Desktop Mechanical Heated Ultrasonic Cleaner Series YR05234



Product Description

Ultrasonic Cleaner Machine YR05234 Series Ultrasonic Cleaner Machine has a various capacity of internal operation, high cleaning accuracy. It can effectively clean the tiny dirt particles, and also clean the deep hole, the small seam and the workpiece concealment.

Product Features

- Proprietary MCU-sweep ultrasonic generator driver circuit to present better cleaning performance and to minimize noise level.
- Proprietary industrial ultrasonic transducer with high Q value and energy conversion efficiency and long service life.
- Imported high-strength glue, staple-free bonding, and long service life of tank.
- Mechanical timed cleaning for 1-60 mins and digital heating control, arbitrary set temperature up to 80 °C.
- Internal tank fabricated of 1.0mm SUS304 stainless steel to present good performance and long service life.
- Cleaning basket made of 304 stainless steel and electro-polished. Outer Shell made of high-quality stainless steel to present good corrosion resistance.



| Model | Internal Dimensions | External Dimensions | Volume (L) | Gross Weight (KG) | Ultrasonic Power(W) | Heating Power(W) | Drain valve | Frequency (KHz) | Time (Min) | Heating Temperature(°C) |
|---------|---------------------|---------------------|------------|-------------------|---------------------|------------------|-------------|-----------------|------------|-------------------------|
| YR05234 | 240*135*100 | 268*170*243 | 3.2 | 4.77 | 120 | 100 | no | | | |
| YR05235 | 240*135*150 | 268*170*290 | 5.0 | 5.4 | 120 | 200 | no | | | |
| YR05236 | 300*150*100 | 325*178*243 | 4.8 | 5.65 | 180 | 200 | no | | | |
| YR05237 | 302*150*150 | 325*178*290 | 6.8 | 6.2 | 180 | 200 | no | | | |
| YR05238 | 300*240*150 | 368*270*290 | 11 | 8.5 | 240 | 300 | yes | 40 | 1-60 | 30-80 |
| YR05239 | 330*300*150 | 397*320*317 | 15 | 11.5 | 360 | 400 | yes | | | |
| YR05240 | 330*300*200 | 397*320*368 | 20 | 12.5 | 360 | 400 | yes | | | |
| YR05241 | 500*300*150 | 560*320*317 | 22 | 14.5 | 480 | 500 | yes | | | |
| YR05242 | 500*300*200 | 560*320*385 | 30 | 16.0 | 600 | 500 | yes | 28 | | |



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW

How does the climate affect your microscope?

microscope is an instrument that allows observing and measuring objects that are too small to be seen with the naked eye. This is achieved through a complex system of lenses and mirrors, which form

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW

What is the importance of the laboratory refrigerator?

refrigerator is a device that is one of the most common pieces of equipment used in laboratories. It is a controlled environment (refrigerated space), various fluids and substances

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89

Cel: +33 (0) 6 80 76 07 10

sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

