




LABORATORY SECTOR

WATHER BATHS CIRCULATOR & CHILLERS



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



KALSTEIN® Wather Baths Circulators & Chillers
Viscometer Circulating Bath
Model YR04985



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



KALSTEIN® Wather Baths Circulators & Chillers
Viscometer Circulating Bath
Model YR04985



Product Description:

YR circulation bath is specially used for viscometers, high temperature control accuracy, low fluctuation.

Features:

- ✓ YR Circulation bath is specially used for viscometers, high temperature control accuracy, low fluctuation
- ✓ Using high quality refrigeration compressor, fast cooling, fully sealed stirring motor, low noise, long serving life.
- ✓ LCD Display can show parameters clearly
- ✓ It has sensor broken protection, prevent malfunction caused by dry burning
- ✓ Special round opening for placing 500 – 250 ml tall beaker
- ✓ The bath also has external circulation, which can direct the controlled liquid to the external thermostatic cup

Model	YR04985	
Specification	Ø75 x 180mm	
Performance	Temperature Range	-5 ~99.9°C
	Temperature Resolution	0.1°C
	Temperature Accuracy	±0.1°C
	Temperature Fluctuation	±0.2°C
Structure	Inner Chamber	High quality Steel plate material
	Outer Case	Corrosion resistant 304 Stainless Steel
	Rated Power	1 000W
	Temperature Control Mode	PID
Controller	Temperature Setting Mode	Touch Button Setting
	Temperature Display	3 Digital led display
	Heating Method	U-shaped heating tuve
	Additional Function	Deviation correction, menú key lock, power off compensation, power off memory, PID adjustable
	Sensor	Thermistor
Safety Device	Over-temperature alarm, short circuit protection, overheat protection, sound-light alarm	
Specification	Inner (WxDxH) (mm)	280 x 140 x 180
	Exterior (WxDxH) (mm)	400 x 400 x 420
	Packing (WxDxH) (mm)	260 x 460 x 450
	Capacity	6 Lts
	Cover Aperture	75mm
	Cover Material	304 Stainless Steel
	Inner Lid Numbere	One
	Net Weight (Kg)	3.2 Kg
Gross Weigh (Kg)	5.6 Kg	



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 

**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:



**Fundación
Jacinto Convit**

<https://www.jacintoconvit.org/>



<https://onetreepanted.org>



Humatem

Solidarité et coopération
biomédicale


<https://www.humatem.org/>



**MANIAPURE
FOUNDATION USA**

<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
Microscopes are instruments used for observing microscopic objects. This is a sensitive device that requires a controlled environment to maintain its performance.
Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
Refrigeration is a key element of the most advanced equipment that is used in laboratories. It is a controlled environment (refrigerated space), without fluctuations in temperature.
Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Alliances
23 views • 3 years ago



Kalstein's Microscopes
445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER
2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

