




LABORATORY LINES

VORTEX MIXER



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



Vortex Mixer

YR05007-1

Description

is best solution for shaking 96 well micro plate and other small size testing tube, centrifuge tube and pcr strip, etc. Applications in cell and molecular biology, microbiology and immunology, etc.

Feature

1. Intelligent controller for speed and timer;
2. Adopt permanent magnet rotor, with over-temperature protection , ensure long time working;
3. Brushless motor, noiseless and maintenance-free;
4. Two operation mode "touching" and "continuous", easy to switch;
5. Can automaticly identify various clamps quickly and ensure safe using;



Input power	AC100-250V/50/60Hz
Output power	18W
Drive mode	DC24V brushless motor
Control mode	PWM steplessly variable frequency
Operation mode	"jog", "continuous"
Power fuse	PPTC Resettable Fuses /no replacement
Maximum load	1.5Kg
Speed Range	0/200-3000r/min.
Oscillation Mode	Circumference
Circle diameter	4.5mm
Fixture Identification function	Automatic identification/ Safe and efficient
Speed function /showing	Full digital / bright LED display
Timing function /showing	Speed timing function /showing Digital timing 1S-999min countdown "jog" timing / bright LED display
Dimension(LxWxH)	198x138x67mm
Weight	4.1kg
The noise for Maximum speed	≤43dB(A)
Maximum acceleration	≤3s
Maximum deceleration	≤3s
Ambient temperature / relative humidity	5-40℃ /80%
Protection grade	IP21
Testing standard	EN 61010-1,EN 61010-2-020,EN 61326-1,EN 61010-3-2/A2



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710/ +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S



Vortex Mixer YR05007-1



Optional working head

 <p>M8-1 standard</p> <p>D < 50mm tubes</p> <p>Max speed: 3000rpm</p>	 <p>M8-3 standard</p> <p>ELISA plate</p> <p>Max speed: 1500rpm</p>	 <p>M8-2 for florence flask</p> <p>D < 85mm flask</p> <p>Max speed: 3000rpm</p>	 <p>M8-4 for testing tube</p> <p>30x10mm</p> <p>Max speed: 1500rpm</p>	 <p>M8-5 for testing tube</p> <p>5x17mm</p> <p>Max speed: 1500rpm</p>	 <p>M8-6 for testing tube</p> <p>16x12-13mm</p> <p>Max speed: 1500rpm</p>
 <p>M8-7 for testing tube</p> <p>12x15mm</p> <p>Max speed: 1500rpm</p>	 <p>M8-8 for testing tube</p> <p>0.2x8ml pcr strip or 96well microplate</p> <p>Max speed: 1500rpm</p>	 <p>M8-9 for testing tube</p> <p>12x17mm</p> <p>Max speed: 1500rpm</p>	 <p>M8-10 for testing tube</p> <p>2x50ml</p> <p>Max speed: 1500rpm</p>	 <p>M8-11 for testing tube</p> <p>4x25mm penicillin bottle</p> <p>Max speed: 1500rpm</p>	



#Letsgivemore

Thanks to your purchase, a donation will be made to a non-profit foundation that fights against breast cancer and helps most vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
Optimum Business Center 450 Rue Baden Powell,
• 34000 Montpellier, France.
Tél: +33 467158849 / +33 680760710 / +33 663810023
<https://kalstein.eu>
KALSTEIN FRANCE, S.A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 


**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:

 **Fundación
Jacinto Convit**
<https://www.jacintoconvit.org/>


 **ONE
TREE
PLANTED**

<https://onetreepanted.org>

 **Humatem**
Solidarité et coopération
biomédicale
<https://www.humatem.org/>

 **MANIAPURE
FOUNDATION USA**
<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed
The 3 articles of the week!
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



1
Read this ARTICLE NOW
How does the climate affect your microscope?
A microscope is an instrument that allows observing and measuring objects that are too small to be seen with the naked eye. This is achieved through a complex system of lenses and mirrors, which form...

Kalstein Blog
Article of the week:
How does the climate affect your microscope?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>



2
Read this ARTICLE NOW
What is the importance of the laboratory refrigerator?
A refrigerator is a device that is one of the most common pieces of equipment used in laboratories to maintain, in a controlled environment (refrigerated space), various fluids and substances...

Kalstein Blog
Article of the week:
What is the importance of the laboratory refrigerator?
Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



MONTPELLIER - FRANCE
HEADQUARTER
Optimum Business Center,
450 Rue Baden Powell,
34000 Montpellier - France.

Fax: +33 (0) 4 67 15 88 49
Tif: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

Official Websites for Latin America.

