




LABORATORY LINES STIRRERS



Kalstein®

A different accompaniment, at your service





Kalstein®

A different accompaniment, at your service

OUR SERVICES

Benefits and Support

In Kalstein France, we take care of the full satisfaction of our customers, that is why we provide value-added services of the highest level based on our experience.



Online Inductions and Trainings

In any part of the world, receive your induction or training from our specialized team of engineers



Quick Response

Our work team is always available to response all your consults or questions, in order to support you in any situation.



#Letsgivemore ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



Technical Support

Enjoy of personalized advice for the correct preventive and corrective maintenance of your equipment, thanks to Kalstein's manuals and articles, special catalogues and video tutorials.



Delivery Logistics

We take care of all the necessary logistics for the dispatch of your goods, whether is by sea, land or air.



Kalstein Worldwide

With more than 25 years growing with our customers, Kalstein's multiformat and modern content, is now present in more than 10 countries and increasing.



Electric Digital Display Overhead Stirrer

Model YR02956




Main Features

- ✓ Brushed DC motor drive , penetration type stirring shaft
- ✓ MCU closed -loop control , constant and stable rotational speed .
- ✓ Digital display and speed setting , rotary encoder operation



Two bases can be chosen:



(U - S H A P E)



(O B L O N G)



#Let'sgive**more** ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



Electric Digital Display Overhead Stirrer Model YR02956



#Letsgivemore ♥

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



Electric Digital Display Overhead Stirrer

Model YR02956



OPTIONAL STIR ROD

Model	Photo	Name	Description
YR400		Stir Rod with shaft	400 mm x Ø 8mm Stainless steel 304 with 3 paddles Diameter 50mm
YR500		Stir Rod with shaft	500 mm x Ø 8mm Stainless steel 304 with 3 -paddles Diameter 50mm
YR400		Stir Rod with shaft PTFE coated	400 mm x Ø 8mm Stainless steel + PTFE
YR500		Stir Rod with shaft PTFE coated	500 mm x Ø 8mm Stainless steel + PTFE .

PACKING LIST

No	Accessory Name	Description
1	Base	320 x 210 x15mm Superfine metal or U - shape base
2	Stir head	170 x 80 x 130 mm. Stainless steel .Motor and Other main Parts
3	Wire	UK plug , EU plug ,etc. 150mm length .
4	Holder	700 mm x Ø 15mm Stainless steel
5	Clamp	200mm length .Can be remove as your requests Stainless steel
6	Stir Rod with shaft	300 mm x Ø 8mm Stainless steel or PTFE material
7	Chuck	50mm length ,stainless steel .



#Letsgivemore

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
2 Rue Jean Lantier • 75001 Paris •
+33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S



Electric Digital Display Overhead Stirrer**Model YR02956**

Model	YR02956
Motor Power	100W.
Motor Type	Brushless DC Motor
Motor Torque	450 m N.M
Machine Input power	120W
Machine output power	100W.
Voltage	220V +/- 10%,or 110V +/- 10%
Rotational Speed Range	60- 1500rpm
Working Time	Continue
Stir Capacity	50ml - 50 L
Stir speed Display	Digital Display
Stir speed Accuracy	+/- 1 rpm .
Stir Rod Length	350mm
Stir Rod Diameter	8mm
Stir Rod Material	Stainless steel 304 .
Holder Length	700mm
Holder Dimeter	25mm
Chuck Clamping Range	Ø1.5-10
Product Dimension (mm)	410 x 93x 180
Gross weight (KG)	12KG
Shell material for machine	Stainless steel 304 with One-time anti-static paint
Viscosity	150, 000mPa.s

**#Let'sgivevmore**

Thanks to your purchase, a donation is made to a non-profit foundation that fights against cancer and helps vulnerable communities.



All rights reserved © KALSTEIN France S. A. S.,
 2 Rue Jean Lantier • 75001 Paris •
 +33 1 78 95 87 89 / +33 6 80 76 07 10 •
<https://kalstein.eu>
KALSTEIN FRANCE, S. A. S





 **Kalstein**[®]
A different accompaniment, at your service

#LETSGIVEMORE 


**WITH THE
ACQUISITION
OF A KALSTEIN
EQUIPMENT**

YOU MAKE A CONTRIBUTION TO:

 **Fundación
Jacinto Convit**
<https://www.jacintoconvit.org/>


 **ONE
TREE
PLANTED**

<https://onetreepanted.org>

 **Humatem**
Solidarité et coopération
biomédicale
<https://www.humatem.org/>

 **MANIAPURE
FOUNDATION USA**
<https://www.maniapurefoundation.org/>

Articles Kalstein
(Social Networks and Official Websites)



Kalstein Blog
Not just be informed, be fully informed

The 3 articles of the week!

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

1
Read this ARTICLE NOW



How does the climate affect your microscope?

microscope is an instrument that allows observing and measuring objects at a very small scale. This is achieved through a complex system of lenses and mirrors, which form

Kalstein Blog

Article of the week:
How does the climate affect your microscope?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

2
Read this ARTICLE NOW



What is the importance of the laboratory refrigerator?

refrigerator is a device that is one of the most common pieces of equipment used in laboratories. It is a controlled environment (refrigerated space), various fluids and substances

Kalstein Blog

Article of the week:
What is the importance of the laboratory refrigerator?

Read more on Kalstein's social medias or official websites
<https://kalstein.eu/articles>

Kalstein Training Videos
(Youtube Channel and Official Websites)



Kalstein France - Jacinto Convit Foundation / Allianc...

23 views • 3 years ago



Kalstein's Microscopes

445 views • 1 week ago



PARIS - FRANCE
HEADQUARTER

2 Rue Jean Lantier
75001 Paris - France

Tlf: +33 (0) 01 78 95 87 89
Cel: +33 (0) 6 80 76 07 10
sales@kalstein.eu

<https://kalstein.eu/>
<https://www.kalstein.fr/>

